


## C30 Hopanes

### Occurrence and origin:

The C30 hopane isomers are normally the most common hopanes of **sedimentary matter**. The origin of the hopanes is the most abundant hopanoid in prokaryotes, C35 tetrahydroxybacteriohopane.



**Cat. No. 0132,30**  
**17 $\alpha$ (H),21 $\beta$ (H)-Hopane**  
 (30 $\alpha\beta$ )



### Geochemical relevance and use in oil spill analysis

Hopanes play an important role in **geochemical investigations**, and are diagnostic biomarker indicators and useful as proof of the origin in **oil spill analysis**, **oil waste analysis** and **analysis of airborne particulates**. They contribute to the so-called terpane fingerprint and are commonly used to relate oils with source rocks.


There are **4 common isomers** of C30-hopanes. The most common are isomers with variable stereochemistry at the 17 and 21 positions, either  $\beta$ (H) with hydrogen above the plane or  $\alpha$ (H) with the hydrogen below the plane.

The **natural isomer  $\beta\beta$**  (17 $\beta$ (H),21 $\beta$ (H)) may be **found in recent sediments**. However, the  **$\alpha\beta$ -isomer** is always the **dominant in mature sediments**, while smaller amounts of the  $\beta\alpha$ - isomer are present. Only minor quantities of the less stable  $\alpha\alpha$ -isomer are present. Thus, the  $\beta\beta$ - and the  $\alpha\alpha$ -isomers are useful internal standards as they normally do not co-elute with other hopanes or triterpenoids in mature sediment.

The  $\beta\alpha$ -isomers (moretanes) are highly specific for **immature to early oil generation**. The moretanes are thermally less stable than the  $\alpha\beta$ -hopanes, and abundances of the C29 and C30 moretanes decrease relatively to the corresponding hopanes with thermal maturity. The ratio of  $\beta\alpha$ -moretanes to their corresponding  $\alpha\beta$ -hopanes decrease with thermal maturity from ca. 0,8 to <0,15. The moretane/hopane ratio is used most commonly for C30, but it is also quantified using C29.

In **fresh oil spills**, the  **$\alpha\beta$ -isomer** of hopane is considered to be non-biodegradable and conserved. Consequently, it can be used as an internal standard to monitor the amount of total oil removed by bioremediation (treatment by oil-degrading bacteria).

The hopanes elute on a normal nonpolar GC-column in the order;  $\alpha\beta$ -,  $\beta\alpha$ -,  $\alpha\alpha$ -,  $\beta\beta$ . The C30 gammacerane (Cat. No. 2646.30) elutes late and in the region **between the C31 22R (1339.31) and C32 22S (1338.31) isomers** while the oleanane isomers ( $\alpha$  and  $\beta$ , Cat. No. **0617.30** and **0618.30**) co-elutes with lupane between  $\beta\alpha$ 29 and  $\alpha\beta$ 30.



*Figure: GC-MS of Mona-2 Oil, Danish North Sea*  
 (Courtesy of Peter Nytoft, GEUS, Denmark)



### Regular C30 Hopanes available from Chiron:

- 5-10 µg neat are supplied in convenient 300µL GC-vials for dilution to e.g. 50-100µg/mL
- 50 and 100 µg/mL are supplied in iso-octane (1 mL ampoules)
- Quantities are measured relative to the intensity (TIC) of 30αβ hopane or by gravimetry

<b>2888.30-50-IO</b>	17α(H),21α(H)-Hopane	50 µg/ml
<b>0132.30-100-IO</b>	17α(H),21β(H)-Hopane	100 µg/ml
<b>0612.30-100-IO</b>	17β(H),21α(H)-Hopane (moretane)	100 µg/ml
<b>0613.30-100-IO</b>	17β(H),21β(H)-Hopane (hopane)	100 µg/ml

#### Other C30 Hopanes

<b>2179.30-50-IO</b>	17a(H),21a (H)-30-Nor-29-methylhopane	50 µg/ml
<b>2262.30-50-IO</b>	17a(H),21β(H)-30-Nor-29-methylhopane	50 µg/ml
<b>2886.30-5UG</b>	17α(H)-30-Diahopane (D30)	5 µg neat
<b>2884.30-5UG</b>	17β(H),21α(H)-22-Methyl-28-nor-spergulane	5 µg neat

#### Other C30 Triterpanes

<b>2646.30-10UG</b>	Gammacerane	10 µg neat
<b>0617.30-100-IO</b>	18α(H)-Oleanane	100 µg/ml
<b>0618.30-100-IO</b>	18β(H)-Oleanane	100 µg/ml
<b>0619.30-100-IO</b>	Friedelane	100 µg/ml
<b>0616.30-100-IO</b>	Lupane	100 µg/ml
<b>0620.30-100-IO</b>	Onocerane I	100 µg/ml
<b>0621.30-100-IO</b>	Onocerane II	100 µg/ml
<b>1192.30-100-IO</b>	20R/20S-Dammarane	100 µg/ml

### Other relevant Biomarker Focus:

**Norhopanes:** Biomarker Focus 7

**Rearranged hopanes:** Biomarker Focus 35

**2-Methyl and 3-Methylhopanes:** Biomarker Focus 37

**Homohopanes and gammacerane:** Biomarker focus 38

#### References:

1. K.E.Peters, C.C. Walters and J.M. Moldowan, *The biomarker guide*, 2. ed. Vol. 1&2, Cambridge University Press, Cambridge 2005.
2. Daling, Faksness, Hansen, and Stout, *Environmental Forensics*, 2002; **3**, 263.
3. cf: <http://www.nordicinnovation.net/nordtestfiler/tec498.pdf>.
4. Wang and Fingas, *Marine Pollution Bulletin*, 2003; **47**, 423, and references therein.
5. Nytoft and Bojesen-Koefoed, *Organic Geochemistry*, 2001; **32**, 841.
6. J.R. Brook *et al.*, *Atmospheric environment*, 2007; **41**, 119-135.

