




Quality Control of Lubricants

In the metal shaping industry lubricant emulsions are used to ease the labor on the piece of work. This is achieved by the emulsion's double function of cooling and lubrication.

The lubricants are divided into two categories:

- non water miscible products
- water miscible products

The water miscible products particularly have to be examined regularly to control the quality of the emulsion determining its performance of cooling and lubrication. The quality of the emulsion also influences the labor process with labor sanitation, health and safety, shelf life of the lubricant and other such costs.




For quality control of the emulsion the following parameters are usually measured:

- | | |
|-------------------------------|------------------|
| * Biozid concentration | * pH |
| * Chloride | * Nitrite |
| * Conductivity | |

At the Daimler - Benz AG plant in Untertürkheim, Stuttgart, engines, axles and gears are manufactured. During this production a variety of different water miscible lubricants are used and analyzed on a regular basis. An important control criterion is the determination of the concentration of the biozid, which is determined by analyzing the parameter Formaldehyde. The standard manual determination of formaldehyde is a very complex and time consuming water distillation procedure whereby about 20 minutes preparation time are required per sample. The weekly number of samples at Daimler - Benz AG amounts to 100, which thus required automation as soon as possible. The Skalar San⁺⁺ Continuous Flow Analyzer provides this automation.

The big advantage of automation in comparison to the classical manual method is the elimination of sample pre-treatment, also the sample volume required is a factor of 3 times less. Specifically for lubricant analysis the random access autosampler can be equipped with a stirring device for homogenizing the sample during uptake.



Lubricant flow in the manufacturing plant



The Skalar San++ Continuous Flow Analyzer



Large range of samplers available

After fine-tuning the Skalar San++ Continuous Flow Analyzer was enabled to run over-night and fully automatically unattended for all installed parameters. At the end of the overnight analysis run the analyzer automatically shuts itself down after completing a rinsing cycle. Results are then available at the appropriate time early in the morning, before production starts. It is then possible to directly correct the concentrations of the essential parameters. An outcome of this was that the reliability of the production process increased significantly.

It should be made clear that parallel to the analysis of the formaldehyde, other important parameters e.g. nitrite, chloride, etc., are also analyzed.

The Skalar San++ Continuous Flow Analyzer has been installed at the Daimler – Benz AG for product quality control and lubricant optimization in the production process. Calibration and concentration correction factors for an optimal lubricant condition were setup in house. These factors fully rationalize the optimum integration effect.

In conclusion the integration of the Skalar San++ Continuous Flow Analyzer has decreased the laboratory labor required for the installed parameters to a factor of 10% of the labor required for the manual methods. The return on investment was calculated at approximately fifteen months on this labor alone. Cost savings on the reduction of reagents used for testing and also in production aspects such as the increment shelf life of the lubricant and less down time during the metal shaping process are also apparent.

Applications for the quality control of lubricants include:

- | | | |
|------------------|----------------|-----------|
| * Aluminium | * Boron | * Calcium |
| * Chloride | * Conductivity | * Cyanide |
| * DOC/TOC | * Formaldehyde | * Iron |
| * Magnesium | * Nitrate | * Nitrite |
| * pH | * Potassium | * Sodium |
| * Total Hardness | * Zinc | |

Skalar's Headquarters
Skalar
 P.O. Box 3237
 4800 DE Breda
 The Netherlands
 Tel. +31 (0)76 5486 486
 Fax. +31 (0)76 5486 400

USA
Skalar, Inc.
 5995 Financial Drive, Suite 180
 Norcross, GA 30071
 Tel. + 1 770 416 6717
 Toll Free: 1 800 782 4994
 Fax. + 1 770 416 6718

Germany
Skalar Analytic GmbH
 Gewerbestraße Süd 63
 41812 Erkelenz
 Germany
 Tel. + 49 (0)2431 96190
 Fax. + 49 (0)2431 961970

France
Skalar Analytique S.A.R.L.
 79, Avenue Aristide Briand
 94110 Arcueil
 France
 Tel. + 33 (0)1 4665 9700
 Fax. + 33 (0)1 4665 9506

United Kingdom
Skalar (UK) Ltd.
 Breda House,
 Millfield Industrial Estate,
 Wheldrake, York, YO19 6NA
 Tel. + 44 (0)1904 444800
 Fax. + 44 (0)1904 444820

Austria
Skalar Analytic GmbH
 Am Anger 22
 A-7451 Oberloisdorf
 Austria
 Tel. + 43 (0)2611 2023411
 Fax. + 43 (0)2611 2023412

Belgium
Skalar Belgium bvba
 Antwerpsestraat 126
 2850 Boom
 Belgium
 Tel. + 32 (0)3888 9672
 Fax. + 32 (0)3844 3441



Email: info@skalar.com
 ©Copyright Skalar 2008

Internet: www.skalar.com
 Publication No. 0106009B. US

Skalar reserves the right to change specifications and appearance of the equipment without further notification.