



## The PRIMACS<sup>SNC</sup> Total Carbon & Total Nitrogen Analyzer




The analysis of Carbon and Nitrogen is important in soil, plant, animal feed, food samples, sediments and sludge samples. The Skalar Primacs<sup>SNC</sup> Analyzer combines the analysis of Total Carbon and Total Nitrogen in one unit using high temperature catalytic combustion.

### Operating Principle

The samples are introduced into the high temperature combustion furnace by the unique vertical “bottom-to-top” sample introduction system. At 1050°C the carbon is completely oxidized to CO<sub>2</sub> in the presence of a catalyst. The CO<sub>2</sub> is measured by Non Dispersive Infra Red Detection (NDIR) for Total Carbon. The analysis of Nitrogen is based upon the well-proven DUMAS technology. Nitrogen is converted in NxOy which is reduced at 600°C to N<sub>2</sub>. The N<sub>2</sub> gas is measured by Thermal Conductivity Detection (TCD). The software displays the carbon and nitrogen peaks simultaneously in real-time and the results can easily be printed or exported to a LIMS system. Whenever priority samples have to be analyzed, the work list can be extended during the run. The sample is weighed into a re-usable quartz crucible and the sample weight is automatically transferred to the work list in the software, which avoids transcription errors.

The Primacs<sup>SNC</sup> analyzer provides an accurate reliable solution for the automation of Total Carbon and Total Nitrogen analysis and has been designed as an easy-to use and low maintenance analyzer. Due to the unique vertical “bottom-to-top” sample introduction, the sample ashes remain in the crucible after the analysis and are taken out of the instrument with removal of the crucible. This avoids sample ash build-up in the combustion zone and therefore reduces the maintenance. The re-usable crucibles and the low consumption of combustion and carrier gas result in a low cost-per-test.

The Primacs<sup>SNC</sup> complies with international regulations such as- CEN, ISO 10694, NEN-EN 131137, AOAC 990.03, AOAC 992.15, AACC 46-30 and ASBC.



## Typical applications:

- Soil and Plant
- Sludges and sediments
- Food
- Waste incinerators
- Waste water plants
- and many more



### General Characteristics

|                     |  |
|---------------------|--|
| Analytes            | Total Carbon and Total Nitrogen / Protein  |
| Method              | TC: High temperature catalytic combustion with Infrared Detection (NDIR)<br>TN/Protein: Dumas combustion technique with Thermal Conductivity Detection (TCD)   |
| Samples             | Solid and liquid samples in - Food, soil, plant, sediment, sludge, waste, etc.   |
| Autosampler         | 20 position integrated autosampler   |
| Sample introduction | Unique vertical "bottom-to-top" sample introduction system   |
| Features            | Automatic balance interfacing<br>Back flush system to remove ambient air for accurate low level analysis<br>Re-usable quartz crucibles<br>Complies with international regulations such as ISO, EN, AOAC, ASBC etc. |

### Operational and Performance Characteristics

|                     |   |
|---------------------|---|
| Measuring range     | Total Carbon: 0.01 – 120 mg C Absolute, Total Nitrogen: 1 – 100 mg N Absolute         |
| Analysis time       | Approx. 3 - 5 minutes   |
| Sample size         | Solid samples: up to 1 gram (100-200 mg nominal)<br>Liquid samples: up to 200 mg (µl) |
| Furnace Temperature | Oxidation furnace 1100°C<br>Reduction furnace 600°C                                   |
| Reproducibility     | < 1 % RSD, FSD (full scale deflection)  |
| Data processing     | Area calculation  |

### Physical Characteristics

|                    |  |
|--------------------|--|
| Gas                | Combustion: O <sub>2</sub> (99.995%), Reduction: He (99.99%) |
| Power requirements | 220 – 240 V, 50/60 Hz  |
| Dimensions (hxdxw) | 83 x 50 x 59 cm (32.7 x 19.7 x 23.2 inches)                  |
| Weight             | 56 kg (123 lb)   |

**Skalar's Headquarters**  
**Skalar**  
 P.O. Box 3237  
 4800 DE Breda  
 The Netherlands  
 Tel. +31 (0)76 5486 486  
 Fax. +31 (0)76 5486 400

**USA**  
**Skalar, Inc.**  
 5995 Financial Drive, Suite 180  
 Norcross, GA 30071  
 Tel. + 1 770 416 6717  
 Toll Free: 1 800 782 4994  
 Fax. + 1 770 416 6718

**United Kingdom**  
**Skalar (UK) Ltd.**  
 Breda House,  
 Millfield Industrial Estate,  
 Wheldrake, York, YO19 6NA  
 Tel. + 44 (0)1904 444800  
 Fax. + 44 (0)1904 444820

**Germany**  
**Skalar Analytic GmbH**  
 Gewerbestraße Süd 63  
 41812 Erkelenz  
 Germany  
 Tel. + 49 (0)2431 96190  
 Fax. + 49 (0)2431 961970

**Austria**  
**Skalar Analytic GmbH**  
 Am Anger 22  
 A-7451 Oberloisdorf  
 Austria  
 Tel. + 43 (0)2611 2023411  
 Fax. + 43 (0)2611 2023412

**France**  
**Skalar Analytique S.A.R.L.**  
 79, Avenue Aristide Briand  
 94110 Arcueil  
 France  
 Tel. + 33 (0)1 4665 9700  
 Fax. + 33 (0)1 4665 9506

**Belgium**  
**Skalar Belgium bvba**  
 Antwerpsestraat 126  
 2850 Boom  
 Belgium  
 Tel. + 32 (0)3888 9672  
 Fax. + 32 (0)3844 3441



Email: [info@skalar.com](mailto:info@skalar.com)  
 ©Copyright Skalar 2008

Internet: [www.skalar.com](http://www.skalar.com)  
 Publication No. 0406010B. US

Skalar reserves the right to change specifications and appearance of the equipment without further notification.