

Quality is more than a word

ESPEC

Compact Ultra Low Temperature Chamber

MC-711•811



Capable of ultra low temperatures as low as -85°C with our unique instrumentation for program or constant operation.

The Compact Ultra Low Temperature Chamber embodies the high performance, reliability, and durability of a full-size chamber.

The line-up is comprised of a total of four models.

Select either the P-instrumentation for programming temperature cycling or the T-instrumentation for constant operation.

Also choose from two very wide temperature ranges that use environmentally-friendly HFC refrigerants.

Select the best model for your specific application and test objectives.





Programming operation type



Constant operation type



Chamber interior



Paperless recorder (optional) *Sample photo

● Four models available

The MC-series comes in two temperature ranges of -75 to $+100^{\circ}\text{C}$ / -85 to $+180^{\circ}\text{C}$, and two types of instrumentation for constant or program operation. A wide temperature range supports tests from temperature characteristic tests to low temperature preservation tests.

● P- and T-instrumentation to meet your test objectives

P-instrumentation with 6.5-inch TFT color LCD enable easy test setting simply by following the displayed instructions. It offers 10 built-in standard programs, and can store up to 20 program patterns (99 steps per pattern), thus capable of diverse program tests. A wide variety of functions such as trend graph display of operation history, timer, and help support are provided for improved operational ease. T-instrumentation with large 7-segment LED offers constant operation.

● Accurate PID temperature control

Just by setting the test temperature, PID control automatically controls temperature, with high accuracy.

● Paperless Recording (optional)

The paperless recorder makes it easy record the temperatures of different components, such as the chamber temperature, on a memory card (Compact Flash).

● Safety measures

Enough precautions are taken to ensure the safety of operators, specimens and the chamber, with various safety measures. In case these safety devices activate, power is shut down to halt chamber operation and details of alarm is displayed on the screen.

● Remote control from your PC

Please contact us for details on using a PC to monitor and remotely control the equipment.

SPECIFICATIONS

Model		MC-711	MC-811
Power supply		200V AC, 3 φ 3W, 50/60Hz	220V AC, 3 φ 3W, 60Hz 380V AC, 3 φ 4W, 50Hz
Maximum current		12A (8.3A at 380V AC)	14A (9.2A at 380V AC)
Temperature control system		Balanced Temperature Control system (BTC system)	
Operating temperature		0 to +40°C (+32 to +104°F)	
Performance ^{*1}	Temperature range ^{*2}	-75 to +100°C (-103 to +212°F)	-85 to +180°C (-121 to +356°F)
	Temperature fluctuation ^{*2}	±0.5°C	
	Temperature uniformity ^{*2}	±1.0°C (±1.8°F)	±1.0°C (±1.8°F): at -85 to +100°C (-121 to +212°F) ±2.0°C (±3.6°F): at +100.1 to +180°C (+212.2 to +356°F)
	Temperature heat-up rate	+20 to +100°C (+68 to +212°F) Approx. 20 min.	+20 to +180°C (+68 to +356°F) Approx. 30 min.
	Temperature pull-down rate	+20 to -70°C (+68 to -94°F) Approx. 60 min.	+20 to -80°C (+68 to -112°F) Approx. 70 min.
Construction ^{*1}	Material	Cold rolled and rust-proof steel plate (Melamine baked finish)	
	Interior	18-8 Cr-Ni stainless steel plate (2B polish)	
	Insulation	Rigid polyurethane foam, fiber reinforced plastics, others	
Heater		Nichrome-stripped wire heater 1kW	
Refrigeration system		Mechanical cascade refrigeration system (air-cooled condenser)	
Refrigerator		Hermetically sealed compressor (R404A/ R508A)	
Refrigerator capacity		650W+400W	800W+650W
Expansion mechanism		Capillary tube system	
Cooler		Plate-fin cooler	
Chamber air circulator		Propeller fan (φ154mm, 4 blades)	
Fittings		Viewing window (φ120mm with frost prevention heater), Cable port (φ50mm, 1pc), Integrating hour meter, Power cable, Drain tube	
Inside dimensions		400W×400H×400D mm (15.7W×15.7H×15.7D in.)	
Outside dimensions ^{*3}		900W×1200H×610D mm (35.4W×47.2H×24.0D in.)	
Capacity		64L (2.2 ft ³)	
Weight		155kg (342 lbs)	

*1 At +23°C (+73.4°F) ambient temperature, no specimen.

*2 The performance is according to JTM K 01-1998 of Japan Testing Machinery Association.

*3 Excluding protrusions.

TEMPERATURE INDICATOR-CONTROLLER

Model		P-instrumentation (SCP-220)	T-instrumentation (ES-102)
Operating mode		Program/ Constant operation	Constant operation
Display		Color TFT LCD display	7-segment LED display
Setting		Analog touch panel method	Mechanical key input
Program memory capacity		RAM pattern: 20 program patterns (99 steps per pattern) ROM pattern: 10 program patterns	—
Setting and indication range	Temp.	-85°C to +110°C (MC-711), -95 to +190°C (MC-811)	
	Time *	0 to 999hrs. 59min.	0 to 99hrs. 59min., 0 to 999hrs.
Setting and indication resolution	Temp.	0.1°C	
	Time *	1 min.	
Input		Thermocouple type T (Copper/ Copper-Nickel)	
Communication function		RS-485	
Auxiliary functions		Input burn-out detection function, Upper and lower temp. limit alarm function, Self-diagnostic function, Alarm indication function, Power cut protection function, Refrigerator capacity automatic control function, Trend graph display function (SCP-220), Help function (SCP-220)	

* For T-instrumentation applies to remote program via the communications functions.

SAFETY DEVICES

- Leakage breaker for power supply (for 200/220V AC only)
- Circuit breaker (for 380V AC only)
- Air circulator temperature switch
- Electric parts compartment door switch
- Control circuit overload & short circuit protection fuse
- Reverse prevention relay
- Refrigerator overload relay
- SSR overload & short circuit protecting circuit breaker
- Thermal fuse
- Specimen power supply control terminal (with power cord plug)
- Compressor temperature switch
- Upper and lower temperature limit alarms (built inside temperature controller)
- Burn-out circuit (built inside temperature controller)
- Watchdog timer (built inside temperature controller)
- Refrigerator automatic delay circuit (built inside temperature controller)
- Overheat protector (independent type)

OPTIONS

- Paperless recorder
- Temperature recorder
- Temperature recorder for future installation
- Thermocouple
- Emergency stop switch *except 380V AC spec.
- External alarm terminal
- Cable port
- Cable port rubber plug
- Shelf / Shelf bracket
- Caster
- Communication function (GPIB/ RS-232C)
- Power cable 5·10m *except 380V AC spec.

ACCESSORIES

- Shelf (stainless)2
- Shelf brackets (stainless) 2 sets
- Cable port rubber plug (φ 50mm)1
- Chamber lamp1
- Glass tube fuse 2 (200, 220 VAC), 1 (380 VAC)
- Thermal fuse1
- Plug type fuse (for 380 VAC)1
- User's manual1



DANGER

● Do not use specimens which are explosive or inflammable, or which contain such substances. To do so could be hazardous, as this may lead to fire or explosion.

● Do not place corrosive materials in the chamber. If corrosive substances or liquid is used, the life of the unit may be significantly shortened specifically because of the corrosion of stainless steel, resin and silicone materials.

● Do not place life forms or substances that exceed allowable heat generation.



CAUTION

Be sure to read the user's manual before operation.

■ Some photographs listed in this catalog contain Japanese display.

ESPEC CORP. <http://www.espec.co.jp/english>

Head Office

3-5-6, Tenjinbashi, Kita-ku, Osaka 530-8550, Japan
Tel: 81-6-6358-4741 Fax: 81-6-6358-5500

ESPEC NORTH AMERICA, INC.

Tel: 1-616-896-6100 Fax: 1-616-896-6150

ESPEC EUROPE GmbH

Tel: 49-89-1893-9630 Fax: 49-89-1893-96379

ESPEC ENVIRONMENTAL EQUIPMENT (SHANGHAI) CO., LTD.

Head Office

Tel: 86-21-51036677 Fax: 86-21-63372237

BEIJING Branch

Tel: 86-10-64627025 Fax: 86-10-64627036

TIANJIN Branch

Tel: 86-22-26210366 Fax: 86-22-26282186

GUANGZHOU Branch

Tel: 86-20-83317826 Fax: 86-20-83317825

SHENZHEN Branch

Tel: 86-755-83674422 Fax: 86-755-83674228

SUZHOU Branch

Tel: 86-512-68028890 Fax: 86-512-68028860

ESPEC TEST TECHNOLOGY (SHANGHAI) CO., LTD.

Tel: 86-21-68798008 Fax: 86-21-68798088

ESPEC (MALAYSIA) SDN. BHD.

Tel: 60-3-8945-1377 Fax: 60-3-8945-1287



ISO 9001/JIS Q 9001
Quality Management System Assessed
and Registered

ESPEC CORP. has been assessed by and registered in the Quality Management System based on the International Standard ISO 9001:2008 (JIS Q 9001:2008) through the Japanese Standards Association (JSA).

ISO 14001 (JIS Q 14001)
Environmental Management System Assessed and Registered
ESPEC CORP.

- Specifications are subject to change without notice due to design improvements.
- Corporate names and trade names mentioned in this catalog are trademarks or registered trademarks.