



A **METTLER TOLEDO** COMPANY

## Prevent Laboratory RSI with Good Technique

## **POOR**







**Seated Posture** 



- Shoulders elevated
- Upper arm elevated
- Elbow extended
- Wrist in deviation



- Lower back supported by chair
- Upper back and neck upright
- Upper arm vertical
- Wrist in the same plane as the forearm

**Standing Posture** 



- Upper back and neck stooped
- Lower back and trunk stooped
- Elbow flexed



- Lower back and trunk upright
- Upper back and neck upright
- Upper arm vertical
- Elbow bent at 90°
- Forearm parallel to the floor
- Wrist in the same plane as the forearm

A2 Good Pipetting Posture.indd 1



- Upper arm flexed
- Elbow extended
- Wrist deviated downward



- Forearm parallel to the floor
- Wrist and forearm in the same plane

17/07/2013 10:28

**Wrist Posture** 



- Wrist extended backwards
- Forearm contact stress on the edge of the bench

Don't Forget these important points when pipetting:-

- Do not over reach
- Work at a comfortable height
- Prevent twisting motions
- Vary pipetting activities
- Alternate with other tasks
- Try using both hands
- Stretch frequently
- Grip the pipette lightly
- Take regular short breaks

For technical papers and other pipetting ergonomics information please visit > www.anachem.co.uk/ergowp