

SiliaMets[®]

Metal Scavengers



Distributed by

Greyhound Chromatography and Allied Chemicals
6 Kelvin Park, Birkenhead, Merseyside CH41 1LT United Kingdom
Tel: +44 (0)151 649 4000 Fax: +44 (0)151 649 4001
sales@greyhoundchrom.com

www.greyhoundchrom.com

Metal Scavenging with SiliaMetS[®]

SiliaBond Metal Scavengers BECOME SiliaMetS Metal Scavengers!

Same Efficient & High Quality Metal Scavenger Products, Brand New Look!

SiliCycle has developed a new look and a new brand for our Metal Scavengers. These products, known as SiliaBond Metal Scavengers (*i.e.*: SiliaBond Thiol), are now named SiliaMetS (*i.e.*, SiliaMetS Thiol) with a new color code. We have updated our branding to give more visibility to our metal removal solutions. This new branding will help differentiate these products from other functionalized silica gels available (*reagents and other bonded phases*).

Although we changed our branding from SiliaBond Metal scavengers to SiliaMetS, no change has been made to the products themselves; you will still be purchasing the same quality products that you have been enjoying for years.

SiliCycle is THE world leader and THE pioneer in metal scavenging solutions. Reasons to choose us:

- Over 12 years of know-how in silica-grafting and metal scavenging technology
- Strong, extensive, and confidential technical support and scientists to help you
- Broadest portfolio of scavengers (*wide variety of ligands*) and applications developed
- Wide range of formats for all purification scales; from laboratory to plant scale purifications
- Cited in many external publications (*and patents*) used by satisfied customers



Introduction

In recent years, the time pressure associated with quickly bringing candidate drugs to market has increased the number of transition metal-catalyzed reactions progressing from lead optimization to early scale-up. The removal of post-reaction metal residues has become a major issue in the pharmaceutical industry. Purification of APIs from residual metal catalyst by traditional methods (*chromatography, activated carbon, distillation, etc.*) often leads to problems such as high costs, time loss, low efficiency, and reduced API yields. To overcome these limitations, SiliCycle has developed **SiliaMetS Metal Scavengers**, a range of products that have significantly changed how chemists purify APIs.

Silica-based metal scavengers have been proven to be the purification method of choice used by

several companies from various industries. Take a look at the section “Customer Success Stories with **SiliaMetS**” to read about our satisfied customers. With the silica matrix advantages over polymers (*no swelling, more general solvent compatibility, higher mechanical and thermal stability, easily scalable applications and availability of different formats, including SPE, flash cartridges, and bulk*) and SiliCycle’s expertise in grafting technology, **SiliaMetS** are the solution of choice for metal removal without contamination of drug candidates. **SiliaMetS** are highly selective and offer a cost-effective alternative for metal removal.


This section includes useful information and tips on **SiliaMetS** (*properties and selection chart*) uses, experimental procedures, and results.



What are SiliaMetS Metal Scavengers?

SiliaMetS Metal Scavengers are functionalized silica gels designed to react and bind excess metal complexes. The process for using scavengers is outlined in the scheme below.

What is a Supported Metal Scavenger?



To be effective, the Metal Scavengers need the ability and inherent functionality to remove metals in their various oxidation states from the reaction mixture. For example, upon completion of a palladium metal-catalyzed reaction, the metal residue contained in the reaction can exist in both Pd (0) and Pd (II).

SiliaMetS - Regulatory Information

For many years, SiliaMetS Metal Scavengers have been used in pilot plants by GMP pharmaceutical, biotechnology, and fine chemical industries as well as contract research and manufacturing organizations. They have run their own analysis proving SiliaMetS Metal Scavengers can safely be used without compromising the purity of the material by leaching of the silica-supported product.

Thus, SiliCycle is committed to high quality standards and always strives to provide defect-free products. In doing so, all products are manufactured in an ISO 9001:2008 compliant facility and subjected to a stringent quality control. Every lot needs to meet the quality specifications and a sample from every batch is kept for subsequent analysis. All products are shipped with the following information:












- Certificate of Analysis
 - Purity (*Leachables and extractables*)
 - Molecular loading
 - Surface Coverage
 - Volatile Content
- Material Safety Data Sheets (*MSDS*)
- BSE/TSE Declaration (*no animal origin*)
- Relevant Technical Information

Need specific regulatory files? SiliCycle can work with you to fill your requirements and provide you custom regulatory documentation including specific analytical tests in line with your needs.

SiliaMetS Product Range

SiliCycle, a leader in functionalized silica gels, has developed a wide range of scavengers to remove a variety of metals at competitive prices.

| SiliaMetS Metal Scavengers Portfolio | | | |
|--------------------------------------|----------------|---|--|
| SiliaMetS | Product Number | Structure | Brief Description |
| SiliaMetS Thiol | R51030B |  | SiliaMetS Thiol is our most versatile and robust metal scavenger for a variety of metals under a wide range of conditions. It has been used in pharmaceutical processes up to production scale. |
| SiliaMetS Thiourea | R69530B |  | SiliaMetS Thiourea is a versatile metal scavenger for all forms of palladium and is widely used in the pharmaceutical industry. Once complexed with a transition metal, it has been reported to be an effective catalyst. |
| SiliaMetS Cysteine | R80530B |  | SiliaMetS Cysteine is the silica bound equivalent of the amino acid cysteine. It is a versatile scavenger for a variety of metals and the preferred metal scavenger for tin residues. By attaching the molecule to the backbone via the amino group, the thiol group remains free and accessible for higher metal scavenging efficiency. |
| SiliaMetS DMT | R79030B |  | SiliaMetS DMT is the silica-bound equivalent of 2,4,6-trimercaptotriazine (<i>trithiocyanuric acid, TMT</i>). It is a versatile metal scavenger for a variety of metals and the preferred metal scavenger for ruthenium catalysts and hindered Pd complexes (<i>i.e. Pd(dppf)Cl2</i>). |
| SiliaBond Amine | R52030B |  | |
| SiliaMetS Diamine | R49030B |  | Better known for their electrophile scavenging efficiency, and their base reagent quality, SiliaMetS Amine, Diamine and Triamine are also proven scavengers for metals. They are very useful for scavenging Pd, Pt, Cr, W and Zn. |
| SiliaMetS Triamine | R48030B |  | |
| SiliaMetS Imidazole | R79230B |  | SiliaMetS Imidazole is a versatile metal scavenger for a variety of metals including, Cd, Co, Cu, Fe, Ni, Pd, Os, and Rh, under a wide range of conditions and the preferred metal scavenger for iron catalysts. |
| SiliaMetS TAAcOH | R69030B |  | SiliaMetS TAAcOH & TAAcONa (<i>Si-Triaminetetraacetic Acid or Sodium Salt</i>) are supported versions of EDTA in their free and sodium salt forms. These two products are effective metal scavengers for Ca, Mg, Li, Ir, Cs, Os, Sn, Pd, Ni and Cu. |
| SiliaMetS TAAcONa | R69230B | | SiliaMetS TAAcOH is effective for metals in low or zero oxidation states, compared to SiliaMetS TAAcONa which is useful for metals in higher oxidation states (<i>2+ or higher</i>). |

All SiliaMetS are made of standard flash silica gel, namely 40 - 63 microns, 60 Å.



| Metals Removed | SilviaMetS Typical Characteristics | | | | | | SilviaMetS |
|---|------------------------------------|------------|-------------------|---------------------|-----------------------------------|--------------------------|----------------------|
| | Color | Endcapping | Molecular Loading | Typical Tap Density | Solvent Compatibility | Prolonged Storage | |
| Ag, Hg, Os, Pd ²⁺ , Pd ⁰ & Ru Cu, Ir, Pd, Rh ⁺ , Rh ²⁺ , Rh ³⁺ , Sc, Sn | White | Yes | 1.20 mmol/g | 682 g/L | All solvents, aqueous and organic | Keep dry | SilviaMetS Thiol |
| Pd ²⁺ , Pd ⁰ Ag, Cu, Fe, Os, Rh ⁺ , Rh ²⁺ , Rh ³⁺ , Sc, Sn | Off-white | Yes | 1.20 mmol/g | 767 g/L | All solvents, aqueous and organic | Keep dry | SilviaMetS Thiourea |
| Cd, Fe, Ir, Os, Ru, Sc & Sn Ca, Cr, Cs, Cu, La, Mg, Pd ²⁺ , Pd ⁰ , Pt, Rh ⁺ , Rh ²⁺ & Zn | Orange | Yes | 0.30 mmol/g | 665 g/L | All organic solvents | Keep dry under argon | SilviaMetS Cysteine |
| Ir, Ni, Os, Pd ²⁺ , Pd ⁰ , Pt, Rh ⁺ , Rh ²⁺ , Rh ³⁺ & Ru Cd, Co, Cu, Fe, Sc & Zn | Light brown | Yes | 0.50 mmol/g | 732 g/L | All solvents, aqueous and organic | Keep dry | SilviaMetS DMT |
| Cd, Cr, Pt, Rh ⁺ & Rh ²⁺ Co, Cu, Fe, Hg, Pb, Pd ²⁺ , W & Zn | Off-white | Yes | 1.20 mmol/g | 700 g/L | All solvents, aqueous and organic | Keep cool (<8°C) and dry | SilviaBond Amine |
| Cr, Pd ²⁺ , Pd ⁰ , & Pt Cd, Co, Cu, Fe, Hg, Ni, Pb, Ru, W & Zn | Off-white | Yes | 1.20 mmol/g | 728 g/L | All solvents, aqueous and organic | Keep cool (<8°C) and dry | SilviaMetS Diamine |
| Cr, Pb, Pd ²⁺ , Pd ⁰ & Pt Co, Cu, Fe, Ni, Ru, W & Zn | Off-white | Yes | 1.20 mmol/g | 736 g/L | All solvents, aqueous and organic | Keep cool (<8°C) and dry | SilviaMetS Triamine |
| Cd, Co, Cu, Fe, Ni, Os, W & Zn Cr, Pd ²⁺ , Pd ⁰ , Rh ⁺ & Rh ²⁺ | Off-white | Yes | 1.20 mmol/g | 681 g/L | All solvents, aqueous and organic | Keep dry | SilviaMetS Imidazole |
| Co, Ni, Os & Sc Cr, Cs, Fe, Pd ²⁺ , Pd ⁰ , Rh ⁺ , Rh ²⁺ & Sn | Off-white | No | 0.40 mmol/g | 635 g/L | All solvents, aqueous and organic | Keep dry | SilviaMetS TAAcOH |
| Cd, Cs, Cu, Fe, Ir, La, Li, Mg, Ni, Os, Rh ³⁺ , Sc, & Sn Cr, Pd ²⁺ , Pd ⁰ , Rh ⁺ , Rh ²⁺ & Zn | Off-white | No | 0.40 mmol/g | 712 g/L | All solvents, aqueous and organic | Keep dry | SilviaMetS TAAcONa |

Preferred SilviaMetS Metal Scavengers for these metals Also Scavenges these metals

Features & Benefits of SiliaMetS Metal Scavengers

SiliaMetS Metal Scavengers are functionalized silica gels designed to react and bind excess metal complexes. The process for using scavengers is outlined in the scheme on page 87.

| Features & Benefits of SiliaMetS | |
|--|---|
| Features | Benefits |
| No leaching | No API contamination by the metal scavenger |
| Very High Purity | Each SiliaMetS product manufactured is submitted to very rigorous quality control in order to provide customers with default-free products and ensure 100% satisfaction |
| High Selectivity | Total recovery of the API |
| Wide Range of Metal Species (<i>various oxidation state</i>) | Efficient for a wide range of metal catalysts |
| Fast Kinetics | Even at room temperature |
| Cost Efficient | Low cost per gram of metal scavenged Less solvent used |
| Solvent Compatibility | Can be used in any solvent, aqueous (<i>pH 2 to 12</i>) and organic |
| New Technologies Compatibility | Suitable for use in microwave synthesizers and flow chemistry |
| Excellent Stability (<i>Thermally and Mechanically</i>) | Works well with overhead stirring Can withstand very high temperatures |
| Ease of Use & Scalable | No swelling or static charge Remove easily by a simple filtration Scalable from mg up to multi-ton scale |
| Various Formats | Amenable to use in SiliaSep & SiliaPrep Cartridges |
| Controlled Loading | Consistent and accurate loading insure lot-to-lot reproducibility |
| Available in Bulk Quantities | Available in large quantities and always in stock |



Metal Scavenging Screening Service

CONFIDENTIALITY
ASSURED

Having a problem removing any residual metal catalyst? Contact us to take advantage of SiliCycle's expertise in metal removal. Our R&D team can find the optimal conditions for you.

Metal Scavenger Screening Services are innovative as they provide an on-hand solution to the pharmaceutical and manufacturing industries. Working with the product that needs to be free of residual metals and the restricted conditions that can be used with the compound (*i.e.*, solvent, temperature), SiliCycle's **Metal Scavenger Screening Service** will quickly develop the most efficient metal scavenging process providing both time and cost savings. Confidentiality is assured, as in most cases the solution involves working with API and other patented materials, and easy technology transfer is guaranteed.

Take the step many major pharmaceutical companies have, and contact us to discuss how we can help you to reach your metal purity goals.

Many screening services adapted to your needs & budgets are available.



SiliaMetS - Typical Experimental Procedures

Screening in Batch Reactor Mode (*bulk*)

To select the best scavenger for initial screening experiments, do the following steps for each SiliaMetS Metal Scavengers included in the kit. Use 4-8 molar equivalents of each SiliaMetS in respect to the residual metal concentration.

1. Dissolve the crude product to be treated in a suitable solvent (*or use directly the crude reaction mixture*) and prepare vials containing the same solution volume.
2. Directly, add each SiliaMetS included in the kit to these vials.
Note: no pre-wetting of the SiliaMetS is required. See “Determining the Optimal Amount of SiliaMetS to use” at page 96.
3. For initial tests, stir the solution for at least one hour at room temperature.
4. Scavenging progress can be followed by normal analytical techniques. The scavenging progress can be estimated by looking at the color of the solution as demonstrated in the figure (*right*). When the scavenging is almost complete, the solution is less colored and SiliaMetS becomes colored. In some occasional cases, if all the samples are still coloured, try one or all of the following: let them react for a longer period of time; add more equivalents of the SiliaMetS, increase the temperature of the reaction.
5. At the end of the scavenging, filter off the SiliaMetS using a fritted funnel or filtration device.
6. Wash the SiliaMetS with additional solvent for total recovery of the API (*or compound of interest*) and concentrate the solution under vacuum.
7. Analyze the residual metal concentration of each vial to identify the most efficient SiliaMetS Metal Scavenger
Note: you can choose more than one scavenger.
8. If you are satisfied with the scavenging efficiency of the best SiliaMetS, direct scale-up is possible. Otherwise, scavenging optimization can be done with SiliaMetS identified in #7 (*see next section*).



Screening with SiliaMetS Fixed Bed Mode (*SPE or Flash Cartridges*)

SiliaMetS fixed bed formats are a great alternative for metal removal and are directly scalable. Initial screening investigations can be done using SiliaPrep 2g/6mL SPE cartridges.

1. Condition the cartridge with 3-4 cartridge volume using the same solvent as the solution to be treated.
2. Add the solution containing the API and the metal to the top of the cartridge and let it pass through the cartridge under gravity.
Note: if needed, a slight positive pressure on the top of the cartridge or a light vacuum can be applied to speed up the flow rate.
3. As shown to the right, a dark coloured band will be observed on the top of the silica bed most of the time.
4. If the residual solution is still coloured, multiple passes through the same cartridge can be done.
5. Once the scavenging is completed, wash the cartridge using at least 3 column volumes of solvent to insure total API (*or compound of interest*) recovery.



Note: in some cases, additional washing may be required.



SiliaMetS Compatible with New Technologies

SiliaMetS In Flow Chemistry

Metal scavenging can also be achieved using SiliaMetS in flow chemistry applications. Simply place SiliaMetS inside the solid-phase reactors provided with your flow system (like Syrris Asia® Solid Phase Chemistry Reactors) and let the solution to be purified flow through these reactors. Multiple reactors can be placed in series and reactors can be heated to obtain optimum scavenging results.



SiliaMetS In Microwave

Metal removal using SiliaMetS can also be done under microwave irradiation to provide excellent scavenging efficiency in just minutes. Simply mix the scavenger and the API dissolved in a suitable solvent a microwave tube and set-up the system with the appropriate parameters. Usually, 5 minutes is enough to scavenge all residual metals.



Experiment Optimization with SiliaMetS

If, upon completion of the screening procedure, the scavenging is not complete or you wish to either reduce the number of equivalents or the reaction time, optimization steps can be undertaken.


Various parameters can be changed one at a time or simultaneously to improve the metal removal efficiency.

Note: you can mix multiple SiliaMetS to get superior efficiency.

Number of SiliaMetS Equivalents

For initial screening experiments we suggest 4-8 molar equivalents be used in respect to the residual metal concentration of each SiliaMetS. Once the preferred scavenger is identified, further optimization can be done to reduce the number of equivalents used (*typically down to 2-4 equivalents*).

Graph represents residual concentration (%) of $\text{Pd}(\text{OAc})_2$ with SiliaMetS Thiol in DMF.



Subsequent Treatments with SiliaMetS


In some cases (*equilibrium process or the presence of multiple species*), multiple treatments with SiliaMetS is suitable instead of a single treatment with a larger amount.

For optimal results, filtration between each treatment can allow for a higher scavenging efficiency.

Reaction Time

In some cases, where increasing the temperature is impossible, longer contact time with the scavenger can allow higher scavenging efficiency.

Conditions: $\text{Pd}(\text{OAc})_2$, THF, SiliaMetS Thiol, RT.





Temperature



In initial screening, we suggest the scavenging experiments be run at room temperature. Usually, metal scavenging is completed after one hour or so. However, when shorter scavenging times are required, higher scavenging rates can be achieved by

increasing the temperature. SiliaMetS can be safely used at elevated temperature without degradation and can be added either at room temperature or directly to a warm solution.

Solvent

SiliaMetS can safely be used in a wide range of organic and aqueous solvents commonly used in laboratory and in process, such as DMF, DMSO, THF, 2-butanone, alcohols, ethers, chlorinated solvent, etc. As demonstrated in the graphs below, the nature of

the solvent does sometimes influence scavenging efficiency. If scavenging or kinetics are too slow, changing solvent or adding a co-solvent should be considered.



SiliaMetS Format (Mode Used)

One advantage of SiliaMetS is their compatibility with various technologies. They can be used in batch, in fixed bed (*SPE or Flash cartridges*), in flow chemistry,

or in microwave. Scavenging efficiency can be improved by changing the mode used.

Mixing Rate

SiliaMetS are mechanically stable and offer excellent scavenging efficiency in batch processes agitated by overhead and magnetic stirrers, as well as orbital shaking under low to moderate agitation rates.

If required, mixing rates can be increased to get better scavenging results. With faster stirring, you improve SiliaMetS dispersion in solution.

pH of the Aqueous Solution

When the scavenging is done in aqueous solutions, it is possible to use SiliaMetS in a pH range of 2 to 12. Depending on the nature of the SiliaMetS, pH can

modify the functional groups present on the scavengers by charging them. Scavenging can be affected (*i.e., amine groups in acidic conditions*).

Determining the Optimal Amount of SiliaMetS

To get an effective metal removal, the amount of SiliaMetS Metal Scavenger used is very important. You can determine how much scavenger will be needed by one of two ways:

- from the residual concentration (*more accurate method*)
- from the amount of metal catalyst used (*when the residual metal concentration is unknown*)

From residual metal concentration (ppm)

Knowing that the palladium (Pd) level in 800 g of material is 500 ppm (*the oxidation state does not affect the calculation*).

Data needed:

- Loading of the scavenger (SiliaMetS Thiol): 1.2 mmol/g
- Metal molecular weight: Ex. Pd = 106.42 g/mol
- Amount of product to be treated : Ex. 800 g
- Residual concentration of metal: Ex. 500 ppm of Pd

1. Determine the amount of palladium to be scavenged

$$\text{Amount of Pd in mg} = \frac{\text{Residual metal concentration} \times \text{Qty of product to be treated}}{1,000}$$

$$\text{Amount of Pd in mg} = \frac{500 \text{ ppm} \times 800 \text{ g of product}}{1,000} = 400 \text{ mg of Pd in 800 g of product}$$

$$\text{Conversion in mmol of Pd} = \frac{\text{Amount of Pd in mg}}{\text{Metal molecular weight}}$$

$$\text{Conversion in mmol of Pd} = \frac{400 \text{ mg}}{106.42 \text{ g/mol}} = 3.76 \text{ mmol of Pd}$$

2. Calculate the amount of scavenger (SiliaMetS Thiol) to use (1 equivalent)

$$\text{Amount of SiliaMetS Thiol to use} = \frac{\text{Number of mmol of metal concentration}}{\text{SiliaMetS Thiol loading}}$$

$$\text{Amount of SiliaMetS Thiol to use} = \frac{3.76 \text{ mmol of Pd}}{1.2 \text{ mmol/g}} = 3.13 \text{ g of SiliaMetS Thiol for 1 eq.}$$

To scavenge 400 mg of palladium, 3.13 g of SiliaMetS Thiol is needed if using only one equivalent. However, it is highly recommend that a minimum of 4 equivalents be used at first. In this case, the amount of SiliaMetS Thiol will be 4 times higher ($4 \times 3.13 \text{ g} = 12.52 \text{ g}$).

Sometimes, the metal residual concentration is unknown. In such a case, the amount (g) of palladium to be scavenged can be replaced by the amount of metal catalyst used for the reaction:

From amount of metal catalyst used

Data needed:

- Amount of metal catalyst used: Ex. 10 g of Pd(PPh₃)₄
- Metal catalyst molecular weight: Pd(PPh₃)₄ = 1,155.56 g/mol

1. Determine the amount of palladium to be scavenged

$$\text{Amount of Pd in mmol} = \frac{\text{Qty of catalyst used for the reaction used} \times 1,000}{\text{Metal catalyst molecular weight}}$$

$$\text{Amount of Pd in mmol} = \frac{10 \text{ g of Pd(PPh}_3)_4 \times 1,000}{1,155.56 \text{ g/mol}} = 8.65 \text{ mmol of Pd (max to be scavenged)}$$

The amount of SiliaMetS Thiol to be used can then be determined as stated above (*see point 2. above*). In this particular case, one equivalent of SiliaMetS Thiol corresponds to 7.20 g.



SiliaMetS Selection Guide

When selecting a metal scavenger, every parameter must be considered: metal catalyst, solvent, residual reagents, by-products, structure of the API (*or molecule of interest*) and temperature. The following table, shown below, helps customers in selecting the most efficient scavenger for a specific metal and application. However, since some parameters may affect the efficiency of the scavenging, we highly

recommend performing a preliminary screening experiment using the SiliaMetS Metal Scavenger Kit.

SiliCycle also offers a confidential [Metal Scavengers Screening Service](#). Contact us to take advantage of our expertise in metal removal. See page 93 to learn more about this service.

SiliaMetS Metal Scavenger Selection Table

| Scavenger | Ag | Ca | Cd | Co | Cr | Cs | Cu | Fe | Hg | Ir | La | Li | Mg | Ni | Os | Pb | Pd (II) | Pd (0) | Pt | Rh (I) | Rh (II) | Rh (III) | Ru (II) | Sc | Sn | W | Zn |
|--------------|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|---------|--------|----|--------|---------|----------|---------|----|----|---|----|
| Si-Thiol | ■ | | | | | | □ | | ■ | □ | | | | | | ■ | □ | ■ | ■ | | □ | □ | □ | | | | □ |
| Si-Thiourea | □ | | | | | | □ | □ | | | | | | | | □ | ■ | ■ | | □ | □ | □ | | | □ | □ | |
| Si-Cysteine | | □ | ■ | | □ | □ | □ | ■ | | ■ | □ | | □ | | | ■ | □ | □ | □ | □ | □ | | ■ | ■ | ■ | | □ |
| Si-DMT | | | □ | □ | | | □ | □ | | ■ | | | | | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ | □ | | □ |
| Si-Amine | | ■ | □ | ■ | | | □ | □ | □ | | | | | | | □ | □ | | ■ | ■ | ■ | | | | | □ | □ |
| Si-Diamine | | □ | □ | ■ | | | □ | □ | □ | | | | | □ | | □ | ■ | ■ | ■ | | | | | □ | | □ | □ |
| Si-Triamine | | | □ | ■ | | | □ | □ | | | | | | □ | | ■ | ■ | ■ | ■ | | | | | □ | | □ | □ |
| Si-Imidazole | | ■ | ■ | □ | | | ■ | ■ | | | | | | ■ | ■ | | □ | □ | | □ | □ | | | | | ■ | ■ |
| Si-TAAcOH | | | | ■ | □ | □ | | □ | | | | | | ■ | ■ | | □ | □ | | □ | □ | | | | ■ | □ | |
| Si-TAAcONa | ■ | ■ | | □ | ■ | ■ | ■ | ■ | | ■ | ■ | ■ | ■ | ■ | ■ | | □ | □ | | □ | □ | ■ | □ | ■ | ■ | | □ |

■ Preferred scavengers □ Scavenges







SiliaMetS Selection Guide (con't)

SiliaMetS Metal Scavengers Selection Guide (Only Catalyst in Solution)

| SiliaMetS | Catalyst, Solvent & Conditions (% of catalyst scavenged) | | | | | | | |
|---------------------|--|--|------------------------------------|------------------------------------|--------------------------|-----------------------------|-----------------------------|--------------------------------|
| | Pd(OAc) ₂ | Pd ₂ (allyl) ₂ Cl ₂ | Pd ₂ (dba) ₃ | Pd(PPh ₃) ₄ | PdCl ₂ (dppf) | Grubbs 1 st Gen. | Grubbs 2 nd Gen. | Hoveyda-Grubbs 1 st |
| | DMF 4 eq., 4 h, 22°C | DMF 4 eq., 4 h, 80°C | DMF 4 eq., 4 h, 22°C | DMF 4 eq., 4 h, 80°C | DMF 4 eq., 4 h, 22°C | DMF 8 eq., 16 h, 80°C | DMF 8 eq., 16 h, 80°C | DMF 8 eq., 16 h, 80°C |
| SiliaMetS Thiol | > 99 | > 99 | 98 | 98 | | 96 | 99 | 93 |
| SiliaMetS Thiourea | > 99 | > 99 | 98 | 91 | | 98 | 96 | 98 |
| SiliaMetS Cysteine | not screened | not screened | not screened | 98 | not screened | not screened | not screened | not screened |
| SiliaMetS DMT | 98 | > 99 [22°C] | > 99 | > 99 | Pd: 94, Fe: 92 | > 99 [4 eq.] | 99 [4 eq.] | 98 [4 eq.] |
| SiliaBond Amine | 98 | > 99 | 97 | | | 97 | | |
| SiliaMetS Diamine | > 99 | > 99 | > 99 | 90 | | 99 | 94 | 98 |
| SiliaMetS Triamine | > 99 | 90 | > 99 | 80 | | 95 | | 95 |
| SiliaMetS Imidazole | not screened | not screened | not screened | not screened | | not screened | not screened | not screened |
| SiliaMetS TAAcOH | 98 | 93 | 97 [80°C] | | | | | |
| SiliaMetS TAAcONa | 97 | | 80 [80°C] | | | | | |

Note: other catalysts results are available on request (metal screened but not shown: calcium, cobalt, cesium, copper, iron, iridium, lanthane, tin, & tungsten. Contact us!)

SiliaMetS Metal Scavengers Selection Guide (Catalysts Scavenging in a Reaction)

| SiliaMetS | Catalyst, Solvent, Condition & Reaction | | | |
|---------------------|--|---|---|--|
| | PdCl ₂ (PPh ₃) ₂ , CuI (in DME) | Pd(OAc) ₂ , P(o-tol) ₃ (in i-PrOH, H ₂ O) | RhCl(PPh ₃) ₃ (in Toluene) | FeCl ₃ ·6H ₂ O |
| |  <p>8 eq., 4 h, 22°C Sonogashira Coupling</p> |  <p>5 eq., 4 h, 40°C Suzuki Coupling</p> |  <p>65 eq., 4 h, 22°C Wilkinson Hydrogenation</p> |  <p>5 eq., 4 h, 22°C Michael Addition</p> |
| SiliaMetS Thiol | Pd: 89, Cu: 29 | 98 | | |
| SiliaMetS Thiourea | Pd: 72, Cu: 80 | 92 | 81 | 82 |
| SiliaMetS Cysteine | | 84 | 88 | > 99 |
| SiliaMetS DMT | Pd: 98, Cu: > 99 | > 99 | 94 | 98 |
| SiliaBond Amine | | 80 | 93 | 98 |
| SiliaMetS Diamine | | 80 | | > 99 |
| SiliaMetS Triamine | | | | 98 |
| SiliaMetS Imidazole | | 88 | 92 | 98 |
| SiliaMetS TAAcOH | | | 81 | 98 |
| SiliaMetS TAAcONa | | | 88 | > 99 |

Scavenging > 99 %





Scavenging 95 - 99 %

Scavenging 90 - 94 %

Scavenging 80 - 89 %




| Catalyst, Solvent & Conditions (% of catalyst scavenged) | | | | | | | | SiliaMetS |
|--|--------------------------|-------------------------|-------------------------|--------------------------------------|----------------------------------|---|---|-----------|
| Hoveyda-Grubbs 2 nd | TPAP | Ni(acac) ₂ | Wilkinson's Cat. | [Rh(OAc) ₂] ₂ | H ₂ PtCl ₆ | Pb(OAc) ₂ ·3H ₂ O | Zn(OAc) ₂ ·2H ₂ O | |
| DMF 8 eq., 16 h, 80°C | DCM 4 eq., 16 h, 22°C | DMF 4 eq., 4 h, 22°C | DMF 4 eq., 4 h, 80°C | DMF 4 eq., 4 h, 80°C | DMF 4 eq., 4 h, 80°C | DMF 4 eq., 4 h, 22°C | DMF 4 eq., 4 h, 22°C | |
| | 96 [4 eq.] | | > 99 [16h] | 97 | 80 [16 h] | 97 | > 99 | Thiol |
| | > 99 | | 99 | 97 | | | 97 [80°C] | Thiourea |
| not screened | not screened | 92 | 88 | not screened | 99 | | > 99 | Cysteine |
| 99 [4 eq.] | > 99 [4 eq.] | 97 | > 99 | > 99 | > 99 | 99 | 94 | DMT |
| | > 99 | | > 99 | > 99 | | | > 99 | Amine |
| 90 | 97 [4 eq.] | 99 | > 99 | > 99 [22°C] | > 99 | 81 | > 99 | Diamine |
| 95 | > 99 | 93 | 97 | 97 [22°C] | 97 | > 99 [80°C] | > 99 | Triamine |
| | not screened | 91 [80°C] | 90 | 97 [22°C] | not screened | | > 99 | Imidazole |
| | > 99 [4 eq.] | > 99 | 97 | 96 [16 h] | | | | TAAcOH |
| | > 99 [4 eq.] | > 99 | 88 | > 99 [16 h] | | 90 | > 99 | TAAcONa |

| Catalyst, Solvent, Condition & Reaction | | | | SiliaMetS |
|--|---|---|---|-----------|
| CuCN (in DMF) | iridium Crabtree's Cat. (in DCM) | LaCl ₃ ·LiCl (in DMF) | PhCH ₂ ZnCl (in THF) | |
|  10 eq., 4 h, 22°C Rosemund von-Braun Cyanation |  4 eq., 4 h, 22°C Alkene Hydrogenation |  1 eq., 4 h, 22°C 1,2-Addition on Ketone |  4 eq., 4 h, 80°C Negishi Coupling | |
| 94 | | | | Thiol |
| > 99 | | | | Thiourea |
| > 99 | 86 | Li: 75, La: > 99 | 91 | Cysteine |
| > 99 | | | 84 | DMT |
| 98 | | | 94 | Amine |
| > 99 | | | 95 | Diamine |
| > 99 | | | 91 | Triamine |
| 95 | | | 94 | Imidazole |
| 80 | | | | TAAcOH |
| > 99 | 80 | Li: 95, La: > 99 | 94 | TAAcONa |

SiliaMetS - A GlaxoSmithKline Case Study¹

A metal scavenging study was performed following the synthesis of a key synthetic intermediate obtained by the Suzuki-Miyaura coupling presented below. Various parameters were investigated including the efficiency of SiliaMetS in different formats, scavenging kinetics, intermediate recovery and purity.



Scavenging Efficiency, Recovery & Purity

Small-Scale Scavenging (*Synthesis Scale ~ 5 g*)

Table below shows the most efficient SiliaMetS Metal Scavenger products for the treatment of the reaction mixture after work-up in both bulk and fixed mode bed (*pre-packed SPE cartridges*).

| SiliaMetS Scavenging Efficiency & Intermediate Recovery Results | | | | |
|---|------------------------------------|------------------|---------------------------|-----------------------|
| SiliaMetS | Batch Reactor Mode (<i>Bulk</i>) | | Fixed Mode (<i>SPE</i>) | Intermediate Recovery |
| | 5 eq., 4 h, 22°C | 5 eq., 4 h, 40°C | 6 mL / 1 g | |
| SiliaMetS Thiol | 95% | > 99% | 98% | > 99% |
| SiliaMetS Thiourea | 83% | 93% | 99% | 98% |
| SiliaMetS Cysteine | 84% | 91% | 97% | > 99% |
| SiliaMetS DMT | 97% | > 99% | > 99% | 98% |
| Initial Pd Concentration: | 179 ppm in MTBE | | 76 ppm in Toluene | - |

Scavenging Conclusion

Addition of only 5 equivalents of SiliaMetS products for 4 hours at the end of the reaction reduces the residual metal concentration to single digit ppm.

Recovery & Purity Conclusion

Palladium was completely removed, while the organic compound was not sequestered by SiliaMetS products. No impurities were released.

¹ *Org. Proc. Res. & Dev.*, 2008, 12, 896



Larger Scale Scavenging (*Synthesis Scale ~ 55 g*)

SiliaMetS Metal Scavengers in pre-packed SiliaSep Flash Cartridges are a great alternative for metal removal at process development scale. These cartridges offer excellent scavenging efficiency as

shown by results in t below. After the first run, almost all the palladium is captured. After three runs, less than 1 ppm remained in solution.

| SiliasSep Scavenging Results | |
|------------------------------|------------|
| Run # | Scavenging |
| 1 | 97% |
| 2 | 99% |
| 3 | > 99% |

Initial Pd Concentration: 700 ppm in AcOEt

Experimental Conditions:

Cartridge Size: 120 g of SiliaMetS Thiol
 Nb. Equivalent of SiliaMetS Thiol: 25 eq.
 Solution Volume: 1 liter
 Flow Rate: 40 mL / min



Metal Scavenging in Flow Chemistry (*Preliminary Results*)

Flow chemistry is a relatively new technique that is gaining in popularity for large scale manufacturing because of the small investment needed to be able to produce large quantities in a short time. SiliaMetS Metal Scavengers can also be used in flow chemistry

to scavenge metals. A crude reaction mixture purified using a Syrris ASIA® Flow Chemistry System is presented below.

SiliaMetS Thiol Scavenging Results in Flow Chemistry

| Flow Rate | Solution Volume | Contact Time with SiliaMetS Thiol | Time Needed to Treat the Solution | Scavenging Results |
|-------------|-----------------|-----------------------------------|-----------------------------------|--------------------|
| 1.50 mL/min | 100 mL | 16 min | 1h10 | 94.0% |
| 1.00 mL/min | 100 mL | 24 min | 1h40 | 94.3% |
| 0.75 mL/min | 50 mL | 32 min | 1h10 | 94.5% |
| 0.50 mL/min | 50 mL | 48 min | 1h40 | 95.0% |

Initial Pd Concentration: 547 ppm in EtOAc

Experimental Conditions:

Scavenger Used:
 SiliaMetS Thiol
 SiliaMetS Nb. Equivalent:
 13.5 eq.
 Reactors: 2 x 12 mL
 Reactors in Series

Total Solution Volume:
 100 mL
 Purification Scale: 12.5 g
 Temperature: 22°C



Variation of Phosphorous Ligand Nature & Scavenging

Even for the same metal, a variation in the scavenging efficiency can be observed depending on the nature of the products present in the solution to be treated. For example, the steric hindrance of a catalyst and the electronic effects of the phosphorous ligands, are factors influencing the removal of the metal. The same Suzuki coupling shown on page 27 was performed using

different phosphorous ligands; three monodentate and three bidentate ligands. For comparison purposes, scavenging screening was done by using the same two sets of conditions. No optimization was done to increase SiliaMetS performance. By experience, using longer reaction times or higher temperatures will allow for better results.

SiliaMetS Scavenging Results with Monodentate Ligands

| SiliaMetS | Triphenylphosphine [PPh ₃] | | Tri(o-tolyl)phosphine [P(otol) ₃] | | Tri-n-butylphosphine [PnBu ₃] | |
|---------------------------|--|------------------|---|------------------|---|------------------|
| | 4 eq., 4 h, 22°C | 4 eq., 4 h, 60°C | 4 eq., 4 h, 22°C | 4 eq., 4 h, 60°C | 4 eq., 4 h, 22°C | 4 eq., 4 h, 60°C |
| SiliaMetS Thiol | 70% | 97% | 87% | 96% | 26% | 85% |
| SiliaMetS Thiourea | 55% | 86% | 54% | 82% | 18% | 41% |
| SiliaMetS Cysteine | 69% | 76% | 77% | 90% | 17% | 44% |
| SiliaMetS DMT | 95% | 97% | 95% | > 99% | 36% | 87% |
| Initial Pd Concentration: | 27 ppm in EtOAc | | 84 ppm in EtOAc | | 90 ppm in EtOAc | |

SiliaMetS Scavenging Results with Bidentate Ligands

| SiliaMetS | 1,1'-bis(diphenylphosphino)ferrocene [dppf] | | 1,3-bis(diphenylphosphino)propane [dppp] | | (+/-) BINAP | |
|---------------------------|---|------------------|--|------------------|------------------|------------------|
| | 4 eq., 4 h, 22°C | 4 eq., 4 h, 60°C | 4 eq., 4 h, 22°C | 4 eq., 4 h, 60°C | 4 eq., 4 h, 22°C | 4 eq., 4 h, 60°C |
| SiliaMetS Thiol | 50% | 69% | 75% | 90% | 31% | 56% |
| SiliaMetS Thiourea | 3% | 23% | 40% | 60% | 33% | 21% |
| SiliaMetS Cysteine | 29% | 36% | 47% | 55% | 19% | 29% |
| SiliaMetS DMT | 14% | 22% | 95% | 98% | 41% | 64% |
| Initial Pd Concentration: | 63 ppm in EtOAc | | 93 ppm in EtOAc | | 16 ppm in EtOAc | |

Scavenging Conclusion

In all cases, SiliaMetS DMT and Thiol remained the better scavengers throughout the study, even though there is a variation in the nature of the ligand.

Ruthenium Scavenging with SiliaMetS

Ruthenium-based catalysts are commonly used in organic synthesis, mainly in olefin metathesis reactions [ROM(P) and RCM]. Grubbs and Hoveyda-Grubbs catalysts are the most popular ruthenium-based complexes in this field of applications. Complete ruthenium removal can be tedious using conventional methods.

SiliaMetS allow the maximal tolerated concentration of the residual ruthenium to be reached. A ruthenium scavenging study was conducted and various parameters were investigated in order to learn more about their influence on the scavengers' robustness as well as to establish the best experimental conditions.

Ruthenium Scavenging Results using SiliaMetS


| SiliaMetS | Grubbs 1 st Gen. | | Grubbs 2 nd Gen. | | Hoveyda-Grubbs 1 st Gen. | | Hoveyda-Grubbs 2 nd Gen. | |
|--------------------|-----------------------------|------------------|-----------------------------|------------------|-------------------------------------|------------------|-------------------------------------|------------------|
| | Toluene ¹ | DMF ² | Toluene ¹ | DMF ² | Toluene ¹ | DMF ² | Toluene ¹ | DMF ² |
| SiliaMetS Thiol | 90% | 96% | - | 99% | 97% | 93% | - | - |
| SiliaMetS Thiourea | - | 98% | - | 96% | 97% | 98% | - | - |
| SiliaMetS DMT | 95% | 99% ² | > 99% | 99% ² | > 99% ² | 98% ² | 98% ² | 99% ² |
| SiliaBond Amine | 95% | 97% | 92% | - | - | - | - | - |
| SiliaMetS Diamine | 99% | 99% | 91% | 94% | > 99% | 98% | - | 90% |
| SiliaMetS Triamine | - | 95% | - | - | 93% | 95% | - | 95% |
| SiliaMetS TAAcOH | 93% | - | - | - | - | - | - | - |
| SiliaMetS TAAcONa | 96% | - | 96% | - | 98% | - | - | - |

Exp. Conditions: ¹8 eq. of SiliaMetS, 16 h, 80°C; ²Only 4 eq. of SiliaMetS. Initial concentration: 500 ppm for all ruthenium-based catalysts.

Note: SiliaMetS Cysteine and Imidazole were not screened in this study (*and are not currently available for this application*). Only SiliaMetS results higher than 90% are presented in this table.

SiliaMetS vs Other Purification Methods

The use of SiliaMetS to remove ruthenium catalyst after a ring-closing metathesis (RCM) reaction is the most effective purification method. As demonstrated below, the main advantage is that no product is lost during the purification step.



Scavenging Results for Various Purification Methods¹

| Scavenging | Scavenger | Filtration over packed bed of ... ² | | | Flash Purification | |
|---------------------|----------------------------|--|--------|--------|--------------------|----------------|
| | SiliaMetS DMT ¹ | Act. Carbon | Celite | Silica | Manual | SiliaSep Cart. |
| Ruthenium captation | 93% | 73% | 24% | 58% | 70% | 73% |


¹ Using 4 eq., 16h, 22°C. ² Solution is passed directly on a packed bed of various adsorbents, which was then washed with the same quantity of solvent.

*Quantitative yield obtained for each purification method (*adjusted in function of the residual concentration of catalyst*). No impurities were generated in all cases using the different methods (*determined by NMR*).

Tin Scavenging with SiliaMetS



Organotin compounds are versatile reagents commonly used in organic synthesis. The two main applications are in Stille couplings or radical reactions. The removal of tin residues can often be an issue due to the high toxicity of this metal.

Traditional removal methods for this impurity are treatment with an aqueous solution of KF, NH₄OH



or NaOH, or with bases such as DBU. However, the efficiency of these methods can vary and may be inapplicable for particular compounds.

SiliaMetS Cysteine & TAAcONa can be used to efficiently remove tin residues from organic mixtures as demonstrated by the examples below.



Tin Scavenging using SiliaMetS Cysteine & TAAcONa

| Reactions | Initial Concentration | SiliaMetS Cysteine | | SiliaMetS TAAcONa | | |
|---------------------------------|-----------------------|---------------------------------|------------------|---------------------------------|------------------|-------------------|
| | | 4 eq., 4 h, 22°C [2 treatments] | 8 eq., 4 h, 22°C | 4 eq., 4 h, 22°C [2 treatments] | 8 eq., 4 h, 22°C | 4 eq., 16 h, 22°C |
| Stille coupling #1 ¹ | 3,385 ppm | 99% | 64% | 96% | 62% | - |
| Stille coupling #2 ¹ | 981 ppm | 90% | 66% | 66% | 50% | - |
| Radical Reduction | 4,090 ppm | 92% | 88% | 90% | 90% | 90% |

¹ Pd residues were completely removed after only one treatment with SiliaMetS Cysteine.


Osmium Scavenging with SiliaMetS

Osmium products are very useful in organic synthesis. One of the most commonly used is osmium tetroxide (OsO₄), which is a very reliable and powerful reagent for the cis-dihydroxylation of alkenes. However, osmium compounds, in particular OsO₄, are highly poisonous, even at low exposure levels, and must be handled with appropriate precautions.

Therefore, it is important to efficiently remove residual osmium from products of interest.

A scavenging study on three organic reactions involving osmium reactants were performed. The metal scavenging efficiency of SiliaMetS is highlighted in the table on the following page.

Osmium Scavenging with SiliaMetS (con't)



Osmium Scavenging using SiliaMetS

| SiliaMetS | Dihydroxylation | | Sharpless Dihydroxylation | | Lemieux-Johnson Oxidation | |
|---------------------|------------------|------------------|---------------------------|------------------|---------------------------|--|
| | 4 eq., 4 h, 22°C | 8 eq., 4 h, 22°C | 8 eq., 16 h, 22°C | 8 eq., 4 h, 22°C | 8 eq., 16 h, 22°C | |
| SiliaMetS Thiol | 87% | > 98% | > 98% | 87% | 92% | |
| SiliaMetS Cysteine | 89% | > 98% | > 98% | 87% | 91% | |
| SiliaMetS DMT | 92% | 97% | > 98% | 87% | 91% | |
| SiliaMetS Imidazole | 87% | > 98% | > 98% | 89% | 91% | |

Initial Os Concentration:

132 ppm in EtOAc


25 ppm in EtOAc

21 ppm in EtOAc

Note: > 98 % of scavenging means < 0.5 ppm of osmium.

Multiple Metal Scavenging with SiliaMetS

SiliaMetS can be used to remove multiple metals in the same reaction with excellent efficiency. The Negishi coupling presented below was performed to show that SiliaMetS can be used to simultaneously remove residual zinc, palladium, and iron present after the reaction.



Multiple Removal Scavenging Results

| SiliaMetS | Palladium | Iron | Zinc |
|---------------------|-----------|-------|-------|
| SiliaMetS Cysteine | 95% | > 99% | 98% |
| SiliaMetS DMT | 83% | 93% | 99% |
| SiliaMetS Imidazole | 84% | 91% | 97% |
| SiliaMetS TAAcONa | 97% | > 99% | > 99% |

Initial Concentration:

188 ppm in THF

110 ppm in THF

6 ppm in THF

Conditions: 4 eq. of SiliaMetS (relative to palladium), 4 h, 22°C.


SiliaMetS Success Stories Published by Customers

SiliaMetS Metal Scavengers are being used by many pharmaceutical companies, several of which are now using them in pilot plants. In the literature, you can find a number of success stories published by customers highlighting the ease of use and reliable performance of **SiliaMetS**. Some examples are presented in the following pages.

An Amgen Case Study¹

In 2009, Amgen published a chapter in "Catalysis of Organic Reactions" related to the use of scavengers for the removal of palladium in small to multi-kilogram production scale. In their study, they evaluated various parameters such as the scavenging efficiency, the influence of the scavenger loading and the loss of product to adsorption (*recovery*). The study was


based on a palladium-catalyzed sulfonamide coupling and scavenger screening was performed at both room temperature and 65°C using 31 different scavengers.




Amgen Scavenger Screening Results

Condition: 20 mg of each scavenger (20% w/w) in 2 mL HPLC vial that contains 1 mL of crude reaction mixture containing 100 mg of product. Each vial was sealed and agitated overnight. Initial palladium concentration was 423 ppm.

The **BEST** scavenger identified during their study was the SiliaMetS Thiourea providing the lowest Pd content (*residual palladium concentration: 3% or < 14 ppm*) without product sequestration. They mentioned that SiliaMetS Thiourea was used extensively in early process development work.



SILIAMETS THIOUREA WAS THE BEST!



SiliaMetS Thiourea

¹Catalysis of Organic Reactions, Chapter 5. Application of Scavengers for the Removal of Palladium in Small Lot Manufacturing


Allgeier & al., Amgen Inc., Thousand Oaks (California)



Cost Comparison for Most Efficient Scavengers ($\geq 80\%$)

At pilot-plant scale, the optimal compromise between the cost per ppm removed and the scavenging efficiency is crucial. The histogram at right shows a cost comparison on best scavengers identified.

Results highlighted by the graph reduced the number of options to only 4 candidates for further evaluation: in pole position the SiliaMetS Thiourea, and then the TMT, TMT-Na₃, and the Smopex 234.



Top 4 Scavengers Overview

A screening validation was conducted on 1 g scale purification (10 mL of solution) with 20% w/w of each top 4 scavengers at 65°C overnight. After filtration, residual metal concentration was analyzed by ICP-MS

and product recovery was determined by HPLC (see below). SiliaMetS Thiourea was chosen for the large scale purification. See Amgen's paper for further details.

Screening Validation Results on Top 4 Scavengers

| Scavengers | Residual Metal Concentration (ppm) | | | Product Recovery | Commentary |
|------------------------|------------------------------------|-----------------------------|---|------------------|---|
| | Screening Exp. in Solution | Validation Exp. in Solution | Validation Exp. in Solid Product ¹ | | |
| SiliaCycle Thiourea | 14 | 11 | 158 | 102% | Best performance but also most expensive. |
| TMT | 33 | 15 | 264 | 104% | Fine in suspension, filterability concerns on scale. |
| Smopex 234 | 36 | 38 | 496 | 84% | Favorable cost but product recovery inadequate |
| TMT-Na ₃ | 85 | 81 | 1 555 | 78% | Very basic compounds (<i>not effective with base-sensitive groups</i>). Low recovery. |
| Purification Scale: | 100 mg | 1 g | 1 g | 1 g | |
| Initial Concentration: | 423 ppm | 381 ppm | 3,577 ppm | | |

Note: ¹Solid product is obtained by dividing the metal concentration in ppm by the amount of product in the test (1g).

Amgen's Conclusion


"Scavengers offer a practical and expedient option for removal of palladium from process streams to ensure quality of organic products... The screening protocol involves treatment of a candidate process stream with 20% w/w scavenger on product at both room temperature and 65°C followed by analysis of Pd and product adsorption. High-temperature treatment increased the efficiency of Pd removal... Evaluation of process costs is key to identifying Pd removal solutions. While scavengers add cost to a process, their use is often justified by the speed to production in early phase development."

An AstraZeneca Case Study

Publication: Ryberg, P., *Organic Process Research & Development*, 12, 2008, 540
Process Chemistry, AstraZeneca PR&D, Sweden.

In 2008, AstraZeneca published a paper in which they removed palladium impurities in a large-scale process. The workup method found to work the best was a treatment with SiliaMetS Thiol (25% w/w or ~1.4 kg)

at 50°C to purify more than 6.7 kg of material. Final residual palladium concentration was as low as 1-2 ppm.




A Pfizer Global R&D Case Study

Publication: Dorow, R.L. & all, *Organic Process Research & Development*, 10, 2006, 493
Pfizer Global Research and Development, Kalamazoo, Michigan (USA)

In 2006, Pfizer published a paper in which they removed palladium & copper impurities in a 20 kg pilot plant scale. They made two subsequent treatments using SiliaMetS Thiol (20% + 7% w/w) at room temperature for 12 hours. After scavenging with SiliaMetS Thiol, the desired product was obtained with a yield of 76% containing only 17 ppm Pd and 1 ppm Cu. An alternative method was also tried using

80% w/w of Deloxan THP (*Degussa AG*) overnight followed by basification with Na₂CO₃. Residual metal concentration with this method was higher compared to that of SiliaMetS and the yield was lower (about 60%-70%). SiliaMetS allows lower residual metal concentration & higher yield with fewer manipulations!






An Abbott Laboratories Case Study

Publication: Ravn, M.M. & all, P., *Organic Process Research & Development*, 14, 2010, 417
Global Pharmaceutical R&D, Process Research & Development and Discovery, Abbott Laboratories, Chicago, Illinois (USA)

In 2010, Abbott Laboratories published a paper in which they removed palladium and iron impurities using SiliaMetS Thiol (50% w/w relative to **1**). Thus,

palladium and iron levels were respectively 6 ppm and 66 ppm. Refer to Abbott's publication for more details.




A Johnson & Johnson Case Study

Publication: Houpis I.N. & all, *Organic Process Research & Development*, 13, 2009, 598
Johnson & Johnson PRD, API Development, Belgium, and Solvias A.G., Synthesis and Catalysis, Switzerland

In 2009, Johnson & Johnson (J&J) in collaboration with Solvias published a paper in which they developed a mild Sonogashira reaction using various metal catalysts. Treatment with SiliaMetS Thiol simultaneously removed Pd, Cu & Al. Residual

concentrations were below 50, 10, and 3 ppm, respectively, in the isolated product **3**. Refer to J&J's publication for more details.

Note: copper comes from a previous synthesis step.




SiliaMetS Leaching & Stability Studies

SiliaMetS Metal Scavengers are being used by many pharmaceutical and biotechnological companies. Each **SiliaMetS** manufactured by SiliCycle is submitted to an extensive washing procedure to ensure the product exhibits extremely low levels of extractables and leachables.

SiliCycle has implemented a quality control procedure to prevent leaching that includes loading and reactivity determination, as well as leachables and extractables analysis (*silica gel purity* $\geq 99.995\%$). The solution must be free of contaminants for the product to successfully pass the rigorous quality control tests.

To address the end users concerns for potential leaching of impurities into reaction mixtures using **SiliaMetS**, we have performed three typical metal containing reactions. We then investigated the detection, identification, and quantification of possible impurities resulting from the scavengers used.

The following three transition metal catalyzed reactions were performed:



Experimental Procedure

Crude reaction mixtures (8 mL) were placed in a standard polypropylene tube equipped with a 20 μ m frit, filled with 1 g of the appropriate **SiliaMetS** Metal Scavenger, and mixed for 4 h at either room temperature or 80°C. Solutions were then filtered through a 0.02 μ m filter prior to analysis.

Leaching Analysis

For each **SiliaMetS**, silane leaching was analyzed by ICP-OES, which has proven to be very sensitive for silicon quantification (*detection limit in solution is 0.125 ppm*). Traces of non-silicon containing impurities were also analyzed by GC-MS and ¹H NMR Analysis. Only results for **SiliaMetS** Thiol and DMT are shown. However, no evidence of impurities was found for all **SiliaMetS**. Contact us for the complete study results.

Gel Purity Calculation Example

$$\text{Impurity \%} = \frac{2 \text{ mg of silicon}}{1,000,000 \text{ mg of SiliaMetS}} \times 100 \Rightarrow 0.0002\% \text{ impurity}$$

$$\text{Gel purity} = 100 - (\text{Impurity \%}) \Rightarrow 99.9998\% \text{ purity}$$



Silane Leaching Analysis by ICP-OES

Results shown in the table below for SiliaMetS Thiol & DMT confirm that minimal leaching occurs with SiliCycle SiliaMetS.


Note: concentration given are in ppm and represent mg of silicon leached per kg of SiliaMetS.

| Stability of SiliaMetS in Suzuki, Stille and Ring-Closing Metathesis reactions | | | | | |
|--|-------------|-----------------|------------|---------------|------------|
| Reaction (solvent) | Temperature | SiliaMetS Thiol | | SiliaMetS DMT | |
| | | [Silicon] | Gel Purity | [Silicon] | Gel Purity |
| Suzuki (Toluene) | 22°C | 2 ppm | 99.9998% | 1 ppm | 99.9999% |
| | 80°C | 2 ppm | 99.9998% | 2 ppm | 99.9998% |
| Stille (1,4-Dioxane) | 22°C | 2 ppm | 99.9998% | 1 ppm | 99.9999% |
| | 80°C | 1 ppm | 99.9999% | 3 ppm | 99.9997% |
| Ring-Closing Met. (DCM) | 22°C | 2 ppm | 99.9998% | 2 ppm | 99.9998% |

Note: Very low levels of silicon were detected in most experiments, giving product purities higher than 99.995%.


Non-Silicon Leaching Analysis

Gas chromatography-mass spectrometry (GC-MS)



Compared to the silica blank spectrum (*bare silica in solvent*), neither experiment showed evidence of any impurities for either SiliaMetS Thiol or DMT.

¹H NMR Analysis (d₆-dms_o)



Note: each experiment was run on a 1 g aliquote of SiliaMetS and was shaken for one hour at room temperature. In GC-MS spectrum, peak at 8.5 minutes is the internal standard (*1-fluoronaphthalene, 100 ppm*). In NMR spectrum, peaks at 2.4 and 3.4 ppm are, respectively, d₆-dms_o and water contained in deuterated solvent.

Stability Study (Shelf Life)

SiliCycle certifies that SiliaMetS Metal Scavengers stored under recommended conditions in an undamaged container are guaranteed to perform for over two years from the manufacturing date without loss of performance (*results at right*).

| SiliaMetS Thiol after Two Years | | |
|---------------------------------|---------------|------------|
| Lot # | QC Date | Scavenging |
| 11577 | January 2008 | > 99.9% |
| | October 2010 | 99.6% |
| 12218 | February 2008 | 99.9% |
| | October 2010 | 99.1% |

Scavenging: 1 000 ppm of Pd(OAc)₂ in DMF.

Conditions: 2 eq. of SiliaMetS Thiol, 1 h, 22°C.

SiliaMetS Metal Scavengers Ordering Information

SiliaMetS Bulk Ordering Information

| Metal Scavenger | Part Number | Metal Scavenger | Part Number |
|--------------------|-------------|---------------------|-------------|
| SiliaMetS Thiol | R51030B | SiliaMetS Diamine | R49030B |
| SiliaMetS Thiourea | R69530B | SiliaMetS Triamine | R48030B |
| SiliaMetS Cysteine | R80530B | SiliaMetS Imidazole | R79230B |
| SiliaMetS DMT | R79030B | SiliaMetS TAAcOH | R69030B |
| SiliaBond Amine | R52030B | SiliaMetS TAAcONa | R69230B |

Formats: 5g, 10g, 25g, 50g, 100g, 250g, 500g, 1kg, 5kg, 10kg, 25kg, ... up to multi-ton scale. Contact us for details.

SiliaSep Metal Scavenger Cartridges Ordering Information (see SiliaSep's section at page 157)

| SiliaSep Type Quantity per box | SiliaSep 4 g 2/box | SiliaSep 12 g 1/box | SiliaSep 25 g 1/box | SiliaSep 40 g 1/box | SiliaSep 80 g 1/box |
|-----------------------------------|-----------------------|------------------------|------------------------|------------------------|------------------------|
| SiliaSep Thiol | FLH-R51030B-ISO04 | FLH-R51030B-ISO12 | FLH-R51030B-ISO25 | FLH-R51030B-ISO40 | FLH-R51030B-ISO80 |
| SiliaSep Thiourea | FLH-R69530B-ISO04 | FLH-R69530B-ISO12 | FLH-R69530B-ISO25 | FLH-R69530B-ISO40 | FLH-R69530B-ISO80 |
| SiliaSep Cysteine | FLH-R80530B-ISO04 | FLH-R80530B-ISO12 | FLH-R80530B-ISO25 | FLH-R80530B-ISO40 | FLH-R80530B-ISO80 |
| SiliaSep DMT | FLH-R79030B-ISO04 | FLH-R79030B-ISO12 | FLH-R79030B-ISO25 | FLH-R79030B-ISO40 | FLH-R79030B-ISO80 |
| SiliaSep Amine | FLH-R52030B-ISO04 | FLH-R52030B-ISO12 | FLH-R52030B-ISO25 | FLH-R52030B-ISO40 | FLH-R52030B-ISO80 |
| SiliaSep Diamine | FLH-R49030B-ISO04 | FLH-R49030B-ISO12 | FLH-R49030B-ISO25 | FLH-R49030B-ISO40 | FLH-R49030B-ISO80 |
| SiliaSep Triamine | FLH-R48030B-ISO04 | FLH-R48030B-ISO12 | FLH-R48030B-ISO25 | FLH-R48030B-ISO40 | FLH-R48030B-ISO80 |
| SiliaSep Imidazole | FLH-R79230B-ISO04 | FLH-R79230B-ISO12 | FLH-R79230B-ISO25 | FLH-R79230B-ISO40 | FLH-R79230B-ISO80 |
| SiliaSep TAAcOH | FLH-R69030B-ISO04 | FLH-R69030B-ISO12 | FLH-R69030B-ISO25 | FLH-R69030B-ISO25 | FLH-R69030B-ISO80 |
| SiliaSep TAAcONa | FLH-R69230B-ISO04 | FLH-R69230B-ISO12 | FLH-R69230B-ISO25 | FLH-R69230B-ISO25 | FLH-R69230B-ISO80 |

SiliaSep Metal Scavenger Cartridges Ordering Information

| SiliaSep Type Quantity per box | SiliaSep 120 g 2/box | SiliaSep 220 g 1/box | SiliaSep 330 g 1/box | SiliaSep XL 800 g 1/box | SiliaSep XL 1600 g 1/box |
|-----------------------------------|-------------------------|-------------------------|-------------------------|----------------------------|-----------------------------|
| SiliaSep Thiol | FLH-R51030B-IS120 | FLH-R51030B-IS220 | FLH-R51030B-IS330 | FLH-R51030B-IS750 | FLH-R51030B-I1500 |
| SiliaSep Thiourea | FLH-R69530B-IS120 | FLH-R69530B-IS220 | FLH-R69530B-IS330 | FLH-R69530B-IS750 | FLH-R69530B-I1500 |
| SiliaSep Cysteine | FLH-R80530B-IS120 | FLH-R80530B-IS220 | FLH-R80530B-IS330 | FLH-R80530B-IS750 | FLH-R80530B-I1500 |
| SiliaSep DMT | FLH-R79030B-IS120 | FLH-R79030B-IS220 | FLH-R79030B-IS330 | FLH-R79030B-IS750 | FLH-R79030B-I1500 |
| SiliaSep Amine | FLH-R52030B-IS120 | FLH-R52030B-IS220 | FLH-R52030B-IS330 | FLH-R52030B-IS750 | FLH-R52030B-I1500 |
| SiliaSep Diamine | FLH-R49030B-IS120 | FLH-R49030B-IS220 | FLH-R49030B-IS330 | FLH-R49030B-IS750 | FLH-R49030B-I1500 |
| SiliaSep Triamine | FLH-R48030B-IS120 | FLH-R48030B-IS220 | FLH-R48030B-IS330 | FLH-R48030B-IS750 | FLH-R48030B-I1500 |
| SiliaSep Imidazole | FLH-R79230B-IS120 | FLH-R79230B-IS220 | FLH-R79230B-IS330 | FLH-R79230B-IS750 | FLH-R79230B-I1500 |
| SiliaSep TAAcOH | FLH-R69030B-IS120 | FLH-R69030B-IS220 | FLH-R69030B-IS330 | FLH-R69030B-IS750 | FLH-R69030B-I1500 |
| SiliaSep TAAcONa | FLH-R69230B-IS120 | FLH-R69230B-IS220 | FLH-R69230B-IS330 | FLH-R69230B-IS750 | FLH-R69230B-I1500 |



SiliaMetS Metal Scavengers Ordering Information (con't)

SiliaSep OT Metal Scavenger Cartridges (rated 60 psi)

| Silica Weight Quantity per box | 2 g 20/box | 5 g 20/box | 10 g 16/box | 15 g 16/box | 20 g 16/box |
|-----------------------------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| SiliaSep OT Thiol | SPE-R51030B-12U | SPE-R51030B-20X | FLH-R51030B-70Y | FLH-R51030B-70i | FLH-R51030B-70Z |
| SiliaSep OT Thiourea | SPE-R69530B-12U | SPE-R69530B-20X | FLH-R69530B-70Y | FLH-R69530B-70i | FLH-R69530B-70Z |
| SiliaSep OT Cysteine | SPE-R80530B-12U | SPE-R80530B-20X | FLH-R80530B-70Y | FLH-R80530B-70i | FLH-R80530B-70Z |
| SiliaSep OT DMT | SPE-R79030B-12U | SPE-R79030B-20X | FLH-R79030B-70Y | FLH-R79030B-70i | FLH-R79030B-70Z |
| SiliaSep OT Amine | SPE-R52030B-12U | SPE-R52030B-20X | FLH-R52030B-70Y | FLH-R52030B-70i | FLH-R52030B-70Z |
| SiliaSep OT Diamine | SPE-R49030B-12U | SPE-R49030B-20X | FLH-R49030B-70Y | FLH-R49030B-70i | FLH-R49030B-70Z |
| SiliaSep OT Triamine | SPE-R48030B-12U | SPE-R48030B-20X | FLH-R48030B-70Y | FLH-R48030B-70i | FLH-R48030B-70Z |
| SiliaSep OT Imidazole | SPE-R79230B-12U | SPE-R79230B-20X | FLH-R79230B-70Y | FLH-R79230B-70i | FLH-R79230B-70Z |
| SiliaSep OT TAAcOH | SPE-R69030B-12U | SPE-R69030B-20X | FLH-R69030B-70Y | FLH-R69030B-70i | FLH-R69030B-70Z |
| SiliaSep OT TAAcONa | SPE-R69230B-12U | SPE-R69230B-20X | FLH-R69230B-70Y | FLH-R69230B-70i | FLH-R69230B-70Z |

SiliaSep OT Metal Scavenger Cartridges (rated 60 psi)

| Silica Weight Quantity per box | 25 g 10/box | 50 g 10/box | 70 g 10/box | 100 g 12/box |
|-----------------------------------|-----------------|-----------------|-----------------|------------------|
| SiliaSep OT Thiol | FLH-R51030B-95K | FLH-R51030B-95M | FLH-R51030B-95N | FLH-R51030B-276F |
| SiliaSep OT Thiourea | FLH-R69530B-95K | FLH-R69530B-95M | FLH-R69530B-95N | FLH-R69530B-276F |
| SiliaSep OT Cysteine | FLH-R80530B-95K | FLH-R80530B-95M | FLH-R80530B-95N | FLH-R80530B-276F |
| SiliaSep OT DMT | FLH-R79030B-95K | FLH-R79030B-95M | FLH-R79030B-95N | FLH-R79030B-276F |
| SiliaSep OT Amine | FLH-R52030B-95K | FLH-R52030B-95M | FLH-R52030B-95N | FLH-R52030B-276F |
| SiliaSep OT Diamine | FLH-R49030B-95K | FLH-R49030B-95M | FLH-R49030B-95N | FLH-R49030B-276F |
| SiliaSep OT Triamine | FLH-R48030B-95K | FLH-R48030B-95M | FLH-R48030B-95N | FLH-R48030B-276F |
| SiliaSep OT Imidazole | FLH-R79230B-95K | FLH-R79230B-95M | FLH-R79230B-95N | FLH-R79230B-276F |
| SiliaSep OT TAAcOH | FLH-R69030B-95K | FLH-R69030B-95M | FLH-R69030B-95N | FLH-R69030B-276F |
| SiliaSep OT TAAcONa | FLH-R69230B-95K | FLH-R69230B-95M | FLH-R69230B-95N | FLH-R69230B-276F |

SiliaPrep Metal Scavenger Cartridges (see SiliaPrep's section at page 173)

| Formats Quantity per box | 200 mg / 3 mL 50 / box | 500 mg / 3 mL 50 / box | 500 mg / 6 mL 50 / box | 1 g / 6 mL 50 / box | 2 g / 6 mL 50 / box |
|-----------------------------|---------------------------|---------------------------|---------------------------|------------------------|------------------------|
| SiliaPrep OT Thiol | SPE-R51030B-03G | SPE-R51030B-03P | SPE-R51030B-06P | SPE-R51030B-06S | SPE-R51030B-06U |
| SiliaPrep OT Thiourea | SPE-R69530B-03G | SPE-R69530B-03P | SPE-R69530B-06P | SPE-R69530B-06S | SPE-R69530B-06U |
| SiliaPrep OT Cysteine | SPE-R80530B-03G | SPE-R80530B-03P | SPE-R80530B-06P | SPE-R80530B-06S | SPE-R80530B-06U |
| SiliaPrep OT DMT | SPE-R79030B-03G | SPE-R79030B-03P | SPE-R79030B-06P | SPE-R79030B-06S | SPE-R79030B-06U |
| SiliaPrep OT Amine | SPE-R52030B-03G | SPE-R52030B-03P | SPE-R52030B-06P | SPE-R52030B-06S | SPE-R52030B-06U |
| SiliaPrep OT Diamine | SPE-R49030B-03G | SPE-R49030B-03P | SPE-R49030B-06P | SPE-R49030B-06S | SPE-R49030B-06U |
| SiliaPrep OT Triamine | SPE-R48030B-03G | SPE-R48030B-03P | SPE-R48030B-06P | SPE-R48030B-06S | SPE-R48030B-06U |
| SiliaPrep OT Imidazole | SPE-R79230B-03G | SPE-R79230B-03P | SPE-R79230B-06P | SPE-R79230B-06S | SPE-R79230B-06U |
| SiliaPrep OT TAAcOH | SPE-R69030B-03G | SPE-R69030B-03P | SPE-R69030B-06P | SPE-R69030B-06S | SPE-R69030B-06U |
| SiliaPrep OT TAAcONa | SPE-R69230B-03G | SPE-R69230B-03P | SPE-R69230B-06P | SPE-R69230B-06S | SPE-R69230B-06U |