


PureKleen[™] Parts Washer







PureKleen Washer Introduction

Developed over the last 30 years, Suncombe PureKeen Parts Washers are available as single door or two door pass through versions. They are high specification units that employ controllable and recipe variable, high power, spray washing, providing a contained, validatable wash. Accessed via a robust stainless steel insulated hinged door with an inflatable seal, the washing chamber is constructed from 316 stainless steel, with radius-corners of sanitary crevice-free fully drainable construction. The washer provides high energy impingement washing together with low energy flushing chosen by a simple recipe selection, to wash and/or dry Parts and other equipment. Built to a generic design, with the ability to custom engineer to suit clients requirements, they achieve excellent washing results whilst using minimal utilities, providing environmentally friendly, low water and energy usage wash.

Standards

Built to hygienic and sanitary standards and available to comply with ASME BPE, GAMP and 21CFR11, the PureKleen Washers are supplied worldwide to the Biotech, Pharmaceutical, Medical, Healthcare, Personal Care and other critical processing industries. Including washing and drying, they incorporate technologies to provide repeatable patterns to ensure total coverage of all surfaces.

Automation System

The Integrated automation is designed to be operator friendly and simple to use whilst providing flexibility and optimisation. Reliable and robust, they have been developed over the last 20 years, they encompass all elements required to provide a controllable, repeatable automatic system. A range of automation levels are available, starting from entry level, through mid level PLC and HMI versions, to advanced validateable systems. All levels provide a repeatable automated cycle.

See Automation Datasheet for full details

PureKleen PartsWasher



Service Compartment

The Service Compartment is located to the left or right hand of the washer and incorporates the liquid and/or air handling equipment. With maintenance access from the side and/or back.

Loading

The PureKleen Parts Washer is supplied with a trolley for to mount the racks allowing loading and unloading away from the washer.


Internal Rack Loaded



Inflatable Door Seal








Internal Finish



Typical Washer Loaded with Double Rack







Designed Racks



PureKleen Washer Typical PFD

Versatile automation system allows you to build your own recipes from pre-commissioned steps in any order. The variables in each step can be individually configured. The system can either be run as a single pass, so that the wash liquid goes directly to drain after washing or can be used in recirculation mode, so that the sump is filled and recirculates the liquid for a pre-set time.

Dosing Pump

Dosing Lance

with Foot Valve

and Strainer

Client Supply

LSL

Pressure

Relief

Valve

TEMPERATURE

TRANSMITTER

 (π)

Booster Pump

Dosing Pump

Reorder

Level

Low

Level

Client Supplied 10 litre

chemical drum

LSL

OPTIONAL CONDUCTIVITY

TRANSMITTER

OPTIONAL FLOW TRANSMITTER

Reorder

Level

Site Supply

Extraction Hood

Basket 2 Location

Basket 1 Location

PureKleen with Fold

Down Door Option

CONDUCTIVITY

TRANSMITTER

H14

HEPA

Filter

(TT)

Air Supply

Filter

Electrical

Heater

Filter

5)

Standard programmes are included for short clean, long clean and self clean/descale. Additional programmes can be individually configured allowing optimisation to ensure repeatability and maximum efficiency.

*See Automation Datasheet for full details

n	Typical Rinses		(
n e	Cold Water Flush	Single Pass or recirculated	
n	Hot Water Flush	Single Pass or recirculated	I
d	Cold Chemical Rinse	Single Pass or recirculated	
	Hot Chemical Rinse	Single Pass or recirculated	F
	Air purge	To remove water from pipework	
lf ~	Air Blow	Drying air pumped into washing	/
d			1
	Gravity Drain	System Drain	(

	Typical Cycle			
	Pre-Rinse	Deliver cold water as a single pass rinse		
_		90		
	Chemical Rinse	Fill sump with hot chemical — recircu-		
		late		
	Inter-rinse Rinse	Fill sump with cold water — recircu-		
		late		
	Final Rinse	Deliver cold water as a single pass		
		rinse		
	Air purge	Removes water from pipework		
	Air Blow	Air pumped into washing chamber		
	Gravity Drain	System Drain		

Typical Manufacturing Standards*

- Sanitary Construction, fully drainable, crevice free.
- 316 stainless steel contact parts, 304 non-contact parts, T.I.G. Welding

Typical Equipment*

- 316 Stainless Steel pipework
- Large Sanitary Centrifugal Pump
- Sanitary Valves, manual and air operated
- Steam heating
- Calibrated Instruments

LEAN Technology

Adopting LEAN principles, the washers were developed to minimise utilities and wash times, whilst ensuring the safety of the operators and the efficiency of the processing. Our automated systems are configured to incorporate LEAN principles including Overall Equipment Effectiveness, Energy Lean and minimise downtime maintenance.


Testing

All functions of the equipment would be fully wet and dry tested and test results would be documented in the 'Pre -Factory Acceptance Test' (FAT) protocol. Following successful completion of this protocol, the client will be invited to the FAT test, where all tests can be repeated or the pre-FAT tests results can be used.

Validation/ Documentation


The lifecycle approach is adopted (DQ, FDS, HDS, SDS, FAT, SAT, IQ & OQ) with validation being key to every stage of the development process, including Factory Acceptance Testing (FAT), SAT and

PureKleen Washer Layout



Front View

Top View



DOOR 2	DEPTH	DOOR 1	ı
Second Door (if fitted)			HEIGHT B
Side View			Loading Height 800mm

PureKleen Washer Chamber Dimensions

Part #	D	Dimensions mm			
	Width	Depth	Height		
PureKleen 38987/760	760	760	760		
PureKleen 39300/860	760	760	860		
PureKleen 39300/1000	760	760	1,000		

PureKleen Washer Dimensions

Part #	Dimensions mm				
	Width	Depth	Height A	Height B	Door 1/2
PureKleen 38987/760	1,500	800	2,200*	760	700
PureKleen 39300/860	2,000	850	2,400*	860	700
PureKleen 39300/1000	2,000	850	2,400*	1,000	850

^{*} reduce height by 450mm if no drying section required

Typical Utilities		
Water	60 litres/min @ 1 bar	
USP Water (optional)	60 litres/min @ 1 bar	
Air	Minimal use @ 7bar	
Electricity	32 kW 400Vac 3ph+n 50hz	
Steam	500 kg/hr @ 3 bar	
Drain	400 litres/min @ 1 bar	

Common Options		
Second Door	Water pre-heater	
Fold Down Door/s	Conductivity Sensors	
Automatic Door/s	Validation Reports	
Steam Heating	21CFR11 Compliance	
Rinse Aid Addition	Drain Cooler	
Ribbon Printer	Remote Control Panel	

Washing Rack Configuration

Two rack loading each with media delivery points to allow washing on one or two levels. Multiple racks and combinations racks are available

- Standard racks are available with full width baskets
- Standard single or two level baskets
- Special racks are available for hoses, parts, filters etc.
- Custom racks are available which are individually modelled to hold custom parts

^{*}See Manufacturing Standards Datasheet for full details

^{*}See Equipment Standards Datasheet for full details