

Solve your problems of transportation, DDA and Evacuation quickly and efficiently with

## Vario-max, the multifunctional Powered Stairclimber

FREE no obligation onsite demonstration







Even on narrow stairs the **Vario-***max* is reliable. Its innovative climbing mechanism can also handle winding stairs. In addition, it can be used on virtually any type of surface, and. e.g. edges of steps are not damaged.

There are hotels and boarding houses without elevators. For people with a walking disability stairs are a great obstacle that must be conquered. Also transporting loads such as suitcases or laundry containers over stairs is a challenge.

This is precisely the point where the **Vario**-max comes in. As a stair climber for persons or loads, on stairs as well as on even ground, the **Vario**-max is your solution for transportation.

Its low weight and easy handling make the **Vario-**max a transportation aid you don't want to do without anymore

Reduce health risks and economize costs. The **Vario-***max* takes over your heavy jobs.



Manual Handling Solutions 58, Paige Close, The Meadows, Watlington Kings Lynn, Norfolk. PE33 0TQ Tel 01553 811977 Fax 01553 811004 sales@manualhandlingsolutions.co.uk www.manualhandlingsolutions.com

## **Vario-***max*

## The multifunctional transportation aid

- Transportation of persons up to 160 kg
- Transportation of goods up to 170 kg
- Can also be used as a wheelchair
- No damage of the load nor the edges of steps
- High manoeuvrability, even on very narrow stairs

- Can be used on virtually all types of stair coverings
- Automatic braking function when climbing downstairs
- Single step mode
- Failsafe tyres
- Includes full belt safety harness



## Important data

| Technical data  | Vorio-max                  |
|---|----------------------------|
| Maximum load (goods)  | 170 kg                     |
| Maximum load (persons)  | 160 kg                     |
| Climbing speed (continually adjustable)                         | 8 - 23 levels              |
| Rechargeable batteries  | 2 x 12 V / 5 Ah            |
| Capacity with one battery charge (depending on person's weight) | 15 - 25 levels             |
| Height of step  | 21 cm   22.5 cm (extended) |
| Weight of battery pack  | 4.3 kg                     |
| Weight of seat  | 8.2 kg                     |
| Weight of power unit  | 22.5 kg                    |
| Total weight  | 35 kg                      |



Quickly, easily, and without tools you can remove the seat from the stair climber.

























Manual Handling Solutions 58, Paige Close, The Meadows, Watlington, Kings Lynn, Norfolk PE33 0TQ Tel 01553 811977 Fax 01553 811004 sales@manualhandlingsolutions.co.uk www.manualhandlingsolutions.com