




Column Connectors Mini Unions

Glass lined bore - GLT Low dead volume

These unions are specifically designed for connecting flexible silica columns and tubing. The glass lined bore ensures an inert connection and precise inside diameter. Two columns can be connected or a broken column joined with negligible loss in column efficiency or inertness.

Each kit contains a special GLT lined union and all ferrules required for installation.

- Specifically designed to connect capillary columns of similar or dissimilar sizes
- Inert and precise glass lined bore Glass Lined Tubing (GLTTm)
- Low dead volume
- Mini Union ferrules included



Mini Unions

For fused silica capillary columns. Glass lined for inertness. Low thermal mass.

Column ID	2nd Column ID	Pack Contents	Part No.
0.05-0.1mm	0.05-0.1mm	1 Union & 5 Ferrules (0.05-0.1mm ID)	103430
0.1-0.25mm	0.1-0.53mm	1 Union & 5 Ferrules (0.1-0.25mm ID)	103431
0.32mm	0.32-0.53mm	1 Union & 5 Ferrules (0.32mm ID)	103432
0.45-0.53mm	0.45-0.53mm	1 Union & 5 Ferrules (0.45-0.53mm ID)	103433
Replacement Mini Union Nuts		5	103400

Ferrules for Mini Unions (Supplied in packs of 10)

Ferrule Style	Column ID	Part Number
100% Graphite	0.32mm	072630
	0.53mm	072629
15% Graphite / 85% Vespel®	0.05-0.1mm	072695
	0.1-0.25mm	072696
	0.32mm	072697
	0.53mm	072698
	1/16" hole	072693
	Blank (No-hole)	072694
SilTite Metal	0.25mm	073470
	0.32mm	073471
	0.53mm	073472
	1/32"	073473

SilTite™ Mini Unions

For fused silica capillary columns. Glass lined for inertness. Low thermal mass.

Column ID	2nd Column ID	Pack Contents	Part No.
0.1-0.25mm	0.1-0.53mm	1 Union, 2 Nuts & 5 Ferrules (0.1-0.25mm ID)	073550
0.32mm	0.32-0.53mm	1 Union, 2 Nuts & 5 Ferrules (0.32mm ID)	073551
0.45-0.53mm	0.45-0.53mm	1 Union, 2 Nuts & 5 Ferrules (0.45-0.53mm ID)	073552
Replacement Sil	Tite Mini Union Nuts'	5	073553

[†]SilTite nuts must be used with SilTite ferrules

