



Restricted to the ground floor of your home?

Regain the freedom of your home by being able to travel between floors smoothly, quietly and safely either seated in our Compact model or in our Harmony wheelchair lift.

The installation of a through floor lift is a far less complicated operation than you might think. Following the initial survey and the subsequent decision to proceed, our experienced team of fitters will carry out all necessary preparatory work and install the lift in less than two days with very little disturbance to your household.

Sometimes space is limited, therefore the Harmony range features an integral ramp and ingenious car design resulting in a wheelchair lift which is only fractionally larger on the outside than the inside.

For people who do not use a wheelchair, the Harmony Compact seated model requires just two thirds of the space of the standard Harmony.

Operation is simple. Having called the lift from the easily accessible wall control, sited to suit user, the car door can be opened by the same control. The low integral ramp allows easy wheelchair access and a push button inside the car closes the door. A push of the control button will signal the car to travel to the other floor, automatically lifting or replacing the aperture ceiling panel en route, where the reverse procedure is followed to exit the car.



Harmony allows you to move between floors smoothly and safely. Can be used with or without your wheelchair. Fold-up seat available.



Harmony Compact at lower level with fixed seat.

