

AC 500

PHARMACY TABLET COUNTER

The AC 500 offers small scale pharmaceutical and veterinary counting operations with an unparalleled flexibility over a wide range of tablet types..



MACHINE FEATURES

- Counts Tablets - Capsules and Caplets
- Digitally controlled count adjustment.
- Accurate dispensing with excellent repeatability.
- Simple operation with fast product change over.
- Tool free removable product contact parts.
- Constructed in compliance with FDA approved materials.
- Optional full IQ/OQ Validation support documentation
- Small machine footprint

YOUR PARTNER IN PACKAGING



TECHNICAL SPECIFICATION

| | |
|----------------------------|---|
| Output: | Based on short runs a 5mm tablet is up to 1000 per min. |
| Electrical Control System: | KBW LCD Display |
| Electrical Requirements: | 220/240V 50Hz 1 Ph or 110/120V 60Hz power supply supplied. To provide 24V |
| Build Standard: | cGMP, CE marked in accordance with Machinery Directives |
| Product Size: | Tablet Ø5mm - Ø15mm Capsule #5 - #00 |
| Container Size: | 200ml |
| Nett Weight: | 10Kg (21 Lbs) |

OPERATION

- Tablets are loaded onto the table by gradually poring them from a container.
- Adjust the Separator to suit each product by turning the Control Knob.
- The Table feeds the tablets in a circular path, dropping in single file through a light beam.
- A photocell records the passage of the tablet incrementing the count on the LCD display until the desired count has been accumulated.
- At the end of the preset count, the table stops waiting the next command.

www.kbw-packaging.com