

CAN-PIERCER for OXYBABY®

Accessory for the easy measurement of oxygen and carbon dioxide in bottles and cans of various sizes in combination with the analyser OXYBABY®.

Benefits

- minimal gas sample volume required
- reliable test results
- safe and simple operation
- adjustable for different sizes of bottles and cans
- robust design, strong probe
- integrated needle-protection
- for packages with metal- or plastic-lids or caps

Note:

Please find further technical information about the measurement and the data administration on the OXYBABY® data sheet.

Preparation of the OXYBABY® for application as a Can-Piercer:

1. Unscrew the probe protection of the OXYBABY® counterclockwise and take it off.
2. Remove the probe of the OXYBABY® and unscrew the filter counterclockwise. Then screw the filter on the probe holder of the Can-Piercer.
3. Insert the OXYBABY® into the holding brackets of the Can-Piercer.
4. Connect one end of the tube (connector with external thread) to the OXYBABY®, and the other end to the filter on the probe holder.
5. For further details please refer to the operating instructions of the OXYBABY®.

ATTENTION!

The analyser is only to be operated with the special filter!



* analyser and bottle not included in package

Type	Can-Piercer for OXYBABY®
Material	steel
Weight	7.5 kg (without OXYBABY®)
Dimensions (HxWxD)	585 x 155 x 235 mm (23.03 x 6.10 x 9.25 inch) packed, without OXYBABY®
Height of bottles or cans	max. 400 mm (15.75 inch)
Approvals	Company certified according to ISO 9001, ISO 14001 and ISO 22000