

21310

QA4 Surgery Trolley System

Manual Functions

The QA4 Surgery system is both patient trolley and operating table, a total solution for modern day treatment centers, offering exceptional stability, the QA4 is ideal for all aspects of Day surgery including Gynaecology, Orthopaedics, General, ENT and even Ophthalmic work.

Features


Height Adjustment	The height adjustment is by a single pedal independently operated from either side of the trolley
Backrest	The backrest is gas assisted and infinitely variable with a range from 0 to 80°
Trendelenburg	The hydraulic damper tilt system is hand operated from the head end of the trolley and allows 20° of head down tilt (option of 12° reverse trendelenburg available on request). Additional safety feature allows trendelenburg even at lowest height.
X-Ray Translucent	X-Ray translucent patient platform allows 75% C-arm access
Head Section	Removable and interchangeable with specialist options – articulates to +25° / -30° for anaesthetic intubation
Leg Section	Removable ultra lightweight, push button release and self locking mechanism
Side Rail	Stainless steel, British standard rail (Euro available on request) and suitable for all types of operating table attachments (arm tables, lateral positioners, stirrups etc)
Cotsides	The cotsides are polyurethane foam moulded which rotate and are removable. push button release and self locking mechanism
Brakes	Braking system can be operated from either side of the trolley and will active all castors simultaneously; twin wheels ensure a secure and stable operating platform
Steering	5 th Wheel is foot activated from the head end of the trolley, highly manoeuvrable trolley easily guided by one person
Transfusion Pole	Polished Stainless steel, height adjustable fitted with the new innovative quick release and safety locking mechanism – side rail mounted
Cylinder Trough	Moulded construction accepts CD & E Size Cylinders (side loading to base)
Mattress	Manufactured from body moulded foam and through its sensitivity to temperature, this gel feel foam actually moulds to the shape of the patients body, thus relieving pressure by distributing the load as evenly as possible. Developed with the flexible stretch properties of Lycra® this outer material is anti micro bacterial water-resistant, vapour permeable and latex free

Cleaning

Moulded construction and state of the art materials, acrylic capped vacuum form base, Teflon© coated body section and anti bacterial components are all designed for ease of cleaning

Optional Extras

A full range of operating table attachments are available – details attached

Product Specifications**Dimensions**

A – Minimum Height	650mm
A – Maximum Height	1010mm
B – Mattress Length	2040mm
C – Maximum Length	2100mm
D – Mattress Width	600mm
E – Side Rail Width	655mm
F – Cotside Width	770mm
G – Brake Width (Brakes off)	835mm
H – Base Length	1050mm
Mattress Depth	75mm
Castor Diameter	150mm

Weights

Trolley Weight	127kg
Fixed Leg Section (Standard)	6kg
Articulating Leg Section	8kg
Trolley Maximum Weight Limit	160kg
Head Section Weight Limit	25kg
Leg Section Weight Limit	50kg

Range of Movement

Headrest Articulation	+25° / -30°
Backrest Articulation	0 to 80°
Trendelenburg Tilt	20°
Reverse Trendelenburg Tilt	12°
Leg Section Articulation (if fitted)	0 to 45°