

Release: HAPR0112 Release date:immediate

New Product Information

Hiden Gas Analysers at ARABLAB 2015 Visit us on Booth 1011

Hiden are exhibiting their latest laboratory gas analysers at **ARABLAB 2015, 23rd - 26th** March, Dubai UAE. Visit us on Booth 1011. Hiden will feature their latest OGA systems for direct real time analysis, guantification and control of gas related processes ranging in pressure from 100 mbar to 50 bar.

The bench-top mass spectrometer systems are coupled to the user process by the integral robust, flexible, guartz-lined heated capillary with sample consumption rates as low as 1 mL/ min and response times as low as 150 milliseconds. Process interface options enable analyses through the very broad pressure regime and are suited to diverse application areas including top-gas analysis, measurement of dissolved species, multi-stream monitoring of gas feed and process exhaust lines (up to 80 individual gas streams), process/thermal reaction studies, respiratory analysis.

System operation, calibration, data acquisition and display are fully automated, with full manual control selectable when required by the researcher. Integral I/O's permit import of external data, such as temperature and weight, for integration with the mass spectral display and output of process detail for process control functions.

For full details Visit us on Booth 1011 at ARABLAB 2015 or contact Hiden Analytical at: info@hiden.co.uk or visit the main website at www.HidenAnalytical.com.



QGA - Atmospheric Gas Analyser



Close-up: QGA Atmospheric Gas Analyser

Hiden Analytical Ltd. 420 Europa Boulevard Warrington WA5 7UN England ■ +44 [0] 1925 445 225 ■ +44 [0] 1925 416 518 ■ info@hiden.co.uk ■www.HidenAnalytical.com